

# Worktop Light for Powered Hummingbird (LCD & Laptop Models) Quick Installation and Maintenance Instructions



HMA-WT-LEDMNT

## Medical Facility Responsibility

Preventive maintenance checks must be performed regularly to maintain the quality and performance of this product. Any parts that may be broken, missing, worn, distorted, or contaminated in any way should not be used and all affected parts should be replaced immediately. Should the necessity of any repair be suspected; please contact your local distributor.



**WARNING:** It is the responsibility of the end user to ensure all aspects of installation are covered by following the full manual and reviewing warnings/notes found at: <http://www.amico.com/products/workstations-wheels>

## Installation Tools



Scissors

## Installation Reference



Accessory LED  
With Mounting Adhesive



Cable Tie  
x3

## LCD Model Installation Instructions

1. Unlock the worktop using the provided key (lock is located on right side of worktop). Open worktop by lifting from the rear.
2. Open the cord access door and remove the cable management cover. **(Figure 1)**
3. Raise the display to its MAXIMUM height on the display post.
4. Clean the desired mounting location on the bottom of the display.
5. Remove the tape backing from the LED accessory and mount to the desired location on the display. For ideal bonding results, apply moderate pressure to the adhesive section (unfold/open LED) for 15 seconds.
6. Run the LED accessory cable through the cord access door. **(Figure 2)**
7. Plug in the connector and tie the LED cable to the CAT5 cable using the provided cable ties. **(Figure 3)**
8. Replace the cord management cover, close the cord access door. **(Figure 4)**
9. Close and lock the worktop.
10. Refer to Enable LED section below.




Figure 1




Figure 2




Figure 3




Figure 4



Figure 5

# Worktop Light for Powered Hummingbird (LCD & Laptop Models) Quick Installation and Maintenance Instructions



HMA-WT-LEDMNT

## Laptop Model Installation Instructions

1. Unlock the worktop using the provided key (the lock is located on the right side of the worktop). Open the worktop by lifting the worktop from the rear.
2. Plug in the connector and tie the **LED CABLE** to the CAT5 cable using the provided **CABLE TIES**. **(Figure 6)**
3. Place the laptop inside the worktop.
4. Clean the desired mounting location on the laptop.
5. Remove the tape backing from the **LED ACCESSORY** and mount to the desired position on the back of the laptop screen **(Figure 7)**. For ideal bonding results, apply moderate pressure on the **LED ACCESSORY** for 15 seconds. Run excess cable beside/under the laptop.
6. Close and lock the worktop.
7. Refer to Enable LED section below.



Figure 6



Figure 7

## Enable LED

### Standard Control Interface Models

1. Test the light functionality by pressing the light bulb button on the Standard Control Interface.
2. Once the **LED ACCESSORY** illuminates, installation is complete.
3. Adjust the **LED ACCESSORY** to desired angle.

### Advanced Control Interface Models

1. Access login screen by pressing the top left corner of the home screen on the Advanced Control Interface. **(Figure 8)**
2. Enter user name and password and select "Log In". Under "Settings", select "Admin Setup".
3. Select the check box for "Top" in the "Enable Light" section and press "Save" **(Figure 9)**. The main admin screen should display three (3) LED options: "Top", "Middle", and "Bottom".
4. Select "Top" to illuminate the LED. **(Figure 10)**
5. Once the LED illuminates, installation is complete.
6. Adjust the LED to desired angle.



Figure 8




Figure 9



Figure 10

# Worktop Light for Powered Hummingbird (LCD & Laptop Models)

## Quick Installation and Maintenance Instructions



HMA-WT-LEDMNT

### Adjustments

The LED light has a hinge to allow for adjustment to lighting angle.

To adjust brightness on a standard interface, hold the light bulb icon until desired lighting level is reached.

### Troubleshooting Guide

Symptom	Possible Cause	Solution
LED mounting tape is not secure	Incorrect installation	Call distributor to replace part
LED installed but does not illuminate	Improper installation or incomplete software administration setup	Refer to Enable LED section

### Maintenance

It is recommended that once every three (3) months, all locks, bolts, and screws are visually inspected for signs of product wear and tear or damage. Check the product for looseness in the mount or orientation, and ensure it is secure.

### Warranty

During the term of warranty: Within the first twelve (12) months from the date of shipment, Amico Accessories will repair or replace any part which is proven to be defective at no cost. After the twelve (12) month period, Amico Accessories will send the parts to the customer free of charge, however shipping and installation will be borne by the customer.

The warranty is valid only when the product has been properly installed according to Amico Accessories specifications, used in a normal manner, and serviced according to factory recommendations. It does not cover failures due to damage which occur in shipments or failures which resulted from accidents, misuse, abuse, neglect, mishandling, alteration, misapplication or damage that may be attributable to Force Majeure.

#### **AMICO ACCESSORIES DOES NOT HONOR VERBAL STATEMENTS CONCERNING THE WARRANTY.**

The distributor and/or dealer are not sanctioned to create verbal warranties about the product described in this agreement. Any statements will not be honored or be made part of the agreement of sale. This document is the final complete and exclusive terms of the agreement.

#### **THIS WARRANTY IS INCLUSIVE AND REPLACES ALL OTHER WARRANTIES.**

Amico Accessories shall not, under any circumstances be liable for incidental or consequential damages including, but not limited to, profit, loss of sales or injuries to person(s) or property.

Correction of non-compliance as noted above will result in completion of all liabilities of Amico Accessories whether based on agreement, neglect or changed materials, designs or specifications without notice.

All claims for warranty must first be approved by Amico Accessories Customer Service Department: (info@amico-accessories.com or 1-877-264-2697). A valid Return Goods Authorization number must be obtained from Amico Accessories prior to commencement of any warranty claim.