

Operating and Maintenance Instructions

Mira LED Procedure Exam Light Floor Stand



Amico

Contents

Introduction	3
Symbols Used in This Manual	3
Intended Use	4
User Profiles	4
Safety Instructions	4
Scope of Delivery	6
Assembling the Floor Stand	7
Install the Locking Rings	10
Uninstall the Locking Rings	10
Installing the Lighthead	12
Operation	15
Disinfecting, Cleaning, and Sterilization	17
Maintenance	18
Setting the Spring Force	18
Notes on the Arm Systems	19
Selection of Fastening Material	19
Cleaning and Disinfection	19
Testing	19
Disassembling for Servicing	20
Spare Parts	21
Troubleshooting	21
Disassembling for Servicing	21
Technical Data of the Lighthead	22
Technical Data of the Arm System	23
Warranty Policy	24

Introduction

Thank you for your purchase of the Mira LED Procedure Exam Light.

This installation manual gives you instructions on how to install the Mira LED Procedure Exam Light for a Floor Stand version of your purchase. For instructions on operation, maintenance and further description of the product please refer to the Operation and Maintenance Manual included.

Please read the following instructions very carefully and follow the safety instructions and requirements.

If you need information for further instructions on operation, maintenance, and/or descriptions of the product please feel free to contact Amico's Customer Service at: alt-customerservice@amico.com





For technical support please contact: alt-techsupport@amico.com





The product shipped to the job site may not be exactly as shown in the installation drawings.

To learn more about Amico's portfolio of healthcare lighting products, visit <https://www.amico.com/lights>

Symbols Used in This Manual

The warning symbols indicate all instructions that are important for safety. Failure to comply with them can lead to injury, damage to the luminaire or the equipment. In combination with the following signal words the warning symbols means:

 DANGER	Can lead to death or serious injury
 WARNING	Can lead to injury
 CAUTION	Can lead to damage to property
	CE conformity mark
	Comply with operating instructions
	Protective conductor
N	Neutral conductor/return conductor for the electrical current
L1	Live conductor

	Storage humidity
	Storage temperature
	Disposal
	Manufacturer
	Read manual
	Do not step
	Do not lean

Trademarks: All trade names mentioned in these instructions are the exclusive property of the relevant manufacturers

Intended Use

The Mira LED Procedure Exam Light is intended to provide illumination for minor procedures, diagnosis, or examinations. It is to be used near the patient to support treatments, which pose no risk to the patient in the event of any interruption caused by light failure. It is meant for continuous operation.

User Profiles

Medical Staff

All persons with medical training who work in the field for which they were trained.

Cleaning Staff

All persons familiar with national and workplace hygiene requirements.

Electricians

All persons with training in electronics and electrical technology who are familiar with the relevant standards and regulations.

Certified Installers

Persons qualified by training, experience and knowledge of the regulations to install and remove the equipment.

Safety Instructions

Please note that certain duties must be carried out by people with the appropriate corresponding qualifications. The instructions are part of the product and must be kept so they are accessible for all subsequent users:

- These instructions are for the Amico Mira LED Procedure Exam Lights, including floor stand equipment.
- Please carefully read these operating instructions before using the equipment. It will help protect you and others from potential injuries that may occur.
- All work on the light (including repairs) must be done by qualified electricians only. The light must be installed by a certified installer.
- The light must not be modified or manipulated. Use in any other manner than with the original parts may result in different specifications and life-threatening dangers.
- Only persons with medical training are permitted to operate the equipment.
- Only trained cleaning personnel are permitted to clean the equipment.
- The light must not be used in an explosive environment. The power supply to the light is a potential ignition source.
- The light must be used in dry rooms only.
- The light must not be left on unattended.
- Keep all flammable equipment and objects away from the Lighthouse and cone during use.

Safety Instructions

- Do not use a damaged light. Defective cables can also be a potential hazard. Do not position cables near heat sources or sharp edges.
- Eye damage. Never look directly at the light cone.
- Replace damaged lenses before using the light again.
- The light must be connected to the mains with a protective earth to prevent electric shock.
- Never apply additional loads to the Lighthouse and spring arm.
- Do not cover the light with a cloth or similar material while operating.
- If problems are encountered that have not been addressed in these operating instructions, please contact your supplier in the interest of your own safety.

Please keep these operating instructions in a secured place near the device. This ensures easy access to safety instructions and important information should the need arise.

- The inspections and maintenance work, as described on page 18, should only be done by trained service technicians. Anyone else is not authorized to perform inspections and maintenance work.
- The adjustment work should only be carried out by a hospital technician or someone with similar qualifications. Anyone else is not authorized to perform adjustment work.
- This equipment has been built accordingly and is operationally safe. Nonetheless, it may still present some danger, especially if not used according to the operating instructions as specified in this manual.
- Modifications are never to be made for any reason. This is to ensure safety for all persons involved. Any modifications made to the system void the warranty.
- The contents of the operating instructions are subjected to change without notice.

Warnings Levels

DANGER

Indication of hazards that can lead to death or serious injury if measures are disregarded.

WARNING

Indication of hazards that can lead to injury if measures are disregarded.

CAUTION

Indication of hazards that can lead to damage to property if measures are disregarded.

EXPLOSION HAZARD

Non-compliance with this symbol may result in damages or injuries due to explosion.

ELECTRIC SHOCK


Non-compliance with this symbol may result in damages or injuries due to electric shock.

Scope of Delivery


Included Hardware




1x Spring Arm




1x Stand Collar



1x Post




1x Roller Base




1x Corner Cap



1x Power Cord




1x #4 Phillips Screw



2x 1/4" Hex Screws



2x M6 Set Screws



1x Shim






2x Locking Ring

Assembling the Floor Stand

WARNING




DANGER OF INJURY FROM UNFOLDING SPRING ARM

1.	<p>Insert the roll collar to the base cover and add provided setscrew using 3 mm hex key (provided).</p>  <p>Insert the set screws from both sides</p>
2.	<p>Bring the wire cable outside.</p>  <p>Bring out the wire cable from here</p>
3	<p>Insert the base cover with the collar to the floor stand post, raise the base cover to a certain height and tight the set screws. Bring the wire cables from the slotted holes on the post.</p>  <p>Bring the wire cable from this hole</p>

Assembling the Floor Stand

WARNING



DANGER OF INJURY FROM UNFOLDING SPRING ARM

4.	<p>Insert the floor stand post into the roll base and secure it using the supplied Thread-Locking Button Head screws, tightening them with the included 4 mm hex key.</p>	 <p>Insert the 1/4-20 Thread-Locking Button Head screws</p>
5.	<p>Connect the red, black, and green wires from the post cable to the power supply. Next, loosen the collar's setscrews to lower the base cover and align its holes with those on the base plate.</p>	 <p>Loosen the set screws to bring the cover down</p> <p>Connect the wires with power supply</p>
6.	<p>Secure the base cover to the base plate by inserting the four provided 8-32 screws into the designated spots shown in the image, and tighten them using a screwdriver.</p>	 <p>Insert four 8-32 size screws provided</p>

Assembling the Floor Stand

WARNING

DANGER OF INJURY FROM UNFOLDING SPRING ARM


7.	<p data-bbox="219 420 673 525">Remove the plastic cover by unscrewing the screw using a Phillips (star) screwdriver.</p>  <p data-bbox="1112 525 1356 556">Detach the plastic cover</p>
8.	<p data-bbox="219 892 641 924">Slide the spring arm onto the post.</p> 

Assembling the Floor Stand

WARNING

DANGER OF INJURY FROM UNFOLDING SPRING ARM

9. Install the shim and two (2) locking rings.
See below for installing the locking rings.



Install the Locking Rings

1. Separate the multiple coils and inset an end into the groove.
2. Press down around the circumference and wind the ring into the groove.



Uninstall the Locking Rings

1. Insert a screwdriver (or similar flat object) into removal notch.
2. Use tool to pry out first end of ring.
3. Grab ring and spiral it out of the groove.




Assembling the Floor Stand

WARNING


DANGER OF INJURY FROM UNFOLDING SPRING ARM

10. Join the arm connector to the connector in the top of the post.


Ensure correct connector alignment when joining connectors. See label on arm for reference.



Remove label after install. Wipe surface with isopropyl, or similar cleaning solution, to remove adhesive residue.



11. Place the end cap on the post and secure it with the #4 Philips screw.






Installing the Lighthead

CAUTION

DAMAGE CAN BE CAUSED BY FALLING LIGHTHEAD




DANGER OF INJURY THROUGH INCORRECT FITTING OF THE LIGHTHEAD

1.	Remove the screw.	
2.	Rotate the sleeve so the notch and slot are aligned.	
3.	Take out the safety segment.	

Installing the Lighthead

CAUTION

DANGER OF INJURY FROM UNFOLDING SPRING ARM




4.	<p>Take out the plug from the arm.</p> 
5.	<p>Join the arm connector to the lighthead connector.</p> 
6.	<p>Insert Lighthead to the arm.</p> 

Installing the Lighthouse

CAUTION

DAMAGE CAN BE CAUSED BY FALLING LIGHTHEAD

DANGER OF INJURY THROUGH INCORRECT FITTING OF THE LIGHTHEAD

7.	Insert the safety segment back.	
8.	Turn back sleeve and tighten screw according to requirements of the friction from the Lighthouse.	
9.	Clip the handle on the Lighthouse.	

Operation


WARNING

DANGER OF DEATH FROM ELECTRIC SHOCK

- Do not plug in a damaged power cable.
- Replace damaged power cable before using light again.
- The mains voltage and frequency must match the data on the name plate.
- Connect to a mains connection with a protective earth only.
- Do not position the floor-stand such that disconnection of power plug is obstructed.

NOTICE:


- Wind the power cable to the cable retainers when it is not being used.



1. Plug in cable (red circle).
2. Install the cable to the main supply.

TRANSPORT:


- Before transporting the Floor Stand, make sure the lighthead is in the moving position.
- Moving position can be achieved by bringing the lighthead to its lower most position as described.



Operation

Mira LED - All Models

Switch on the Light from the control pad at Lighthouse.



Disinfecting, Cleaning, and Sterilization

NOTICE:

Damage due to improper cleaning.

- Do not use alcohol, solvent or cleaning agents containing chlorine or scouring agents.
- Damage to lens by concentrated disinfectant or solvent.
- Follow the manufacturers recommended mixing ratio.
- Use of incorrect cloths may scratch lens.

NOTICE:

A dirty lens will reduce the illuminating power.

- Clean lens regularly.
- Spray and wipe disinfection is permitted.
- Clean the lens with a leather cloth soaked in cleaner.
- Spray disinfectant in horizontal position.

NOTICE:

The external components of the light can be disinfected with the following diluted products:

- Lysoformin
- Dismozon
- Hexaquart plus
- Sagrotan – fast disinfecting cleaner
- Tuffie5

NOTICE:

The local government regulations for disinfection must also be observed in addition to these instructions.

NOTICE:

Sterilize handle at 120°C and 1.05 bar for 30 minutes.

Maintenance

DANGER

ELECTRIC SHOCKS ARE LIFE-THREATENING


- Remove the connector from the power supply and put the switch in the “OFF” position.
- The connector cable must be checked for damage at least once per year.

NOTES:

- Maintenance and repairs must only be carried out by qualified electricians.
- The corresponding user profile is described in Safety Instructions pages 4-5.

Setting the Spring Force

NOTE: The spring force is set at the optimum setting in the factory.



1. Uninstall the Elbow Caps by removing the screw securing it.
2. Slide the Elbow Cap cover to reveal the adjustment screw and push arm to upper position.
3. Rotate the adjustment screw accordingly to set the spring force:
 - Turn screw counter clockwise (+) to increase spring force (when spring arm lowers).
 - Turn screw clockwise (-) to decrease spring force (when spring arm rises).

Notes on the Arm Systems

Selection of Fastening Material

We recommend the use of (not included in the delivery) counter-plates for fastening to drywall.

Cleaning and Disinfection

DANGER

ELECTRIC SHOCKS ARE LIFE-THREATENING

- The devices can conduct electricity and must be handled with care during cleaning and disinfection.
- If there is a power plug remove it from the socket.
- Do not use spraying methods for cleaning and/or disinfection.
- Do not spray liquid into the openings of electrical power sockets, gas sockets or devices, or allow liquid to enter these.

Recommended Cleaning Agents

Use mild soap solutions or standard detergents as cleaning agents. Wipe down the surfaces of the equipment with a dampened cloth. If necessary, use small quantities of mild soap solution as a cleaning agent. Then wipe the outer surfaces with a soft, clean cloth until they are completely dry.

Disinfection Procedure

Disinfection by wiping is intended to be the standard disinfection procedure for the support arm system. Hygiene guidelines and corresponding safety measures must be determined by the operator for the disinfection procedures that are applied.

- After the contamination by any potentially infectious material (such as blood, secretions or excrement), the areas must be disinfected immediately in a targeted way.
- Comply with application concentrations.
- Wipe the area to be disinfected rather than spraying it.
- Do not use areas that have been wiped until the disinfectant is dry.

Testing

DANGER

ELECTRIC SHOCKS ARE LIFE-THREATENING

During all testing work, disconnect the device from the power supply/remove the plug from the power socket and secure the device so that it cannot be switched on again.

Every 6 months or 5,000 Arm Movements, Check for:

- Paint damage
- Cracks in plastic parts
- Deformation of the support system
- Loosening of parts

Every 2 Years:

- The securing element of the cantilever arm must be checked and greased by a qualified professional.


Notes on the Arm Systems

Disassembling for Servicing



INJURY HAZARD


- The cantilever arm is under a high spring load. If the end device is not moved to the upper cantilever arm position, the cantilever arm will spring up, possibly causing serious injury. Only remove the end device if the cantilever arm is in the upper position.
- Disconnect the device from the power supply/remove the plug from the mains and ensure it cannot be switched on again.
- Place cantilever arm in the upper position.
- Unscrew the brake screw (1) so that the sleeve (2) can be turned.
- Turn the sleeve (2) through 180° until the securing segment (3) appears in the slot (4).
- Hold the end device head (5) firmly and remove the securing segment (3) using a narrow flat-bladed screwdriver.
- Remove end device (5) from the connection.
- Put protective cap in place.




1. Brake screw
2. Sleeve
3. Securing segment
4. Slot
5. End device terminal

This system is designed to be equipped with end devices from other manufacturers. For end device operation instructions please refer to the operating instructions of the related manufacturer.

Spare Parts



Slow-blow 2A Fuse
(L-X-MIRA-FLRSTND-FUSE)



Handle
(L-X-MILED-HDL-2)



Power Cord: US PLUG
(L-X-MILED- FLS-PWR-CRD)

Troubleshooting

Malfunction	Possible Cause	Remedy	User Profiles
Light does not come on	Contact fault	Switch on again	All
	Illuminant defective	Contact manufacturer's After-Sales Service	Only by manufacturer's After-Sales Service
	No mains power	Check power supply and all connections	Qualified electrician
	2A fuse defective (Mira LED Floor Stand only)	Replace fuse	Qualified electrician

Disassembling for Servicing

DANGER

ELECTRIC SHOCKS ARE LIFE-THREATENING

- Remove connector from the power supply and put the switch in the "OFF" position before disassembling.

DANGER

NEVER UNSCREW THE UPPER CABLE HOLDER - INJURY HAZARD

- When both cable holders are unscrewed, the connecting piece is loose and falls downward, which could cause injury and/or damage to the cable and/or the device.



Do not dispose of this product with normal refuse. Dispose of the product at a recycling center or return it to a dealer with a disposal service depending on local regulations. Cut off the cable flush with the housing.

Technical Data of the Lighthouse

	Units	Mira 50 LED	Mira 65 LED	Mira 90 LED
Lighthouse Diameter	in	18	18	18
Light Intensity (Central Luminance)	lx (fc)	50,000 Lux (4,762 fc) @ 1m (39.37")	65,000 Lux (6,190 fc) @ 1m (39.37")	90,000 Lux (8,571 fc) @ 1m (39.37")
Brightness Control	%	50-100	50-100	50-100
Light Field Diameter	in	6.5-12.0	6.5-12.0	6.5-12.0
Depth of Illumination	in	94.5	94.5	94.5
Color Rendering Index	Ra	95	95	95
R9	in	90	90	90
Color Temperature (Standard)	K	4500	4500	4500
Shadow Dilution				
Single Mask	%	0.7	0.7	0.7
Double Mask	%	55.5	55.5	55.5
Cavity	%	100	100	100
Single Mask with Cavity	%	0.3	0.3	0.3
Double Mask with Cavity	%	55.3	55.3	55.3
Number of LED's (Main)		36	36	36
LED Service Life	hours	70,000	70,000	70,000
Total Radiant Power at Maximum Intensity	W/m ²	178	231.4	317
Electrical		100-240 VAC, 50/60 Hz, 0.26 amps	100-240 VAC, 50/60 Hz, 0.31 amps	100-240 VAC, 50/60 Hz, 0.42 amps
Power Supply	ea	Direct Power Supply	Direct Power Supply	Direct Power Supply
Mounting Options		Ceiling, Wall, Floor Stand	Ceiling, Wall, Floor Stand	Ceiling, Wall, Floor Stand
Certificates		cETLus conform to UL60601-1; UL60601-2-41; UL certified IP54 per IEC 60598	cETLus conform to UL60601-1; UL60601-2-41; UL certified IP54 per IEC 60598	cETLus conform to UL60601-1; UL60601-2-41; UL certified IP54 per IEC 60598

Technical Data of the Arm System

Electrical	
Nominal Input Voltage	100-240 VAC
Frequency Range	50 / 60Hz
Input Current	4A-2A
Device Protection	Class 1 Approval ETL/cETL
Any power supply with output less than or equal to 60V, current output less than or equal to 5A and a minimum of two (2) MOOP can be used with the spring arm.	
Environmental Conditions for Transport, Storage, and Operation	
Storage and Transport After 15 weeks the conditions for operation apply.	-10°C to +50°C Ambient 20% to 90% humidity 700 hPa to 1060 hPa
Operation This product should be acclimatized to its operating conditions for 24 hours prior to its first use.	+10°C to +40°C 30% to 75% humidity
Weight	
Mira LED Ceiling Mount With Lighthead	21.5 lbs (9.8 kg)
Operating Mode	
Operating Mode	continuous operation
Classification	
Mira LED Swing	Class I
Classification under 93/42 EEC – Annex IX (medical device class)	Class I
Safety Testing According to:	CSA C22.2#60601-1 IEC/EN 60601-1 2nd. UL 60601-1

Warranty Policy

Medical Lighting Solutions

(Mira Series, Nova Series, Vega Series)

Amico Lights Corporation warrants its products against defective materials and workmanship for twelve (12) months from the date of shipment. Within this period, Amico Lights Corporation will repair or replace any part on site or at the factory which is proven to be defective at Amico Lights Corporation's cost.

Furthermore, Amico Lights Corporation will warrant its material to be free from defect for an additional period of four (4) years (five [5] years from the date of shipment). Within this period, Amico Lights Corporation will replace any part at no charge, which is proven to be defective. Shipping and installation costs after the first twelve (12) months will be borne by the customer.

The spring arms for the above mentioned are warranted for a period of twelve (12) months. The warranty for the spring arms shall include parts and labor and shall commence from the date of shipment.

The warranty applies to normal usage and does not apply to any product that has been subject to alteration, abuse, negligence, or use (including voltage and/or current) other than that for which the product was designed.

Amico Lights Corporation shall not be liable for incidental or consequential damages resulting from the use of the equipment.

All claims for warranty must first be approved by Amico Lights Corporation. A valid Return Goods Authorization (RGA) number must be obtained from Amico Lights Corporation prior to commencement of any service work. Warranty work which has not been pre-authorized by Amico Lights Corporation will not be reimbursed.

Amico Lights Corporation | 122-B East Beaver Creek Road, Richmond Hill, ON L4B 1G6, Canada
Toll Free Phone: 1.877.462.6426 | Toll Free Fax: 1.866.440.4986 | Tel: 905.764.0800 | Fax: 905.764.0862
Email: alt-sales@amico.com | www.amico.com