




DUPLEX SCROLL AIR PLANT WITH DESICCANT AIR TREATMENT SYSTEM (2.0 HP - 7.5 HP)



INLET AND DISCHARGE SIZE VERIFIED

- DIMENSIONS SUBJECT TO CHANGE WITHOUT NOTICE
- ADDITIONAL MODELS AVAILABLE UPON REQUEST

MODEL	HP (kW)	INLET	OUTLET	TANK SIZE (GALLON)	CFM @ 120PSIG	CFM @ 50PSIG	NFPA SYSTEM CAPACITIES WITH ONE RESERVE COMPRESSOR ON STANDBY (@120/@50PSIG)	WEIGHT (LBS)	SQ. FT. REQUIRED (M ²)	SYSTEM BTU/HR	SYSTEM FLA (WITH ALL COMPRESSORS RUNNING)(208V/230V/460V)	SOUND LEVEL dB (A)	LENGTH (IN)	WIDTH (IN)	HEIGHT (IN)
A-SCD-D-080P-SS-N-020	2 (1.5)	1 NPT	0.5 NPT	80	5.6	6.2	5.6 / 6.2	1730 LBS	34.4 (3.20)	5,089	17A/ 15A/ 7A	74	67.00 IN	74.00 IN	83.00 IN
A-SCD-D-080P-SS-N-030	3 (2.2)	1 NPT	0.5 NPT	80	8.5	9.2	8.5 / 9.2	1770 LBS	34.4 (3.20)	7,633	23A/ 20A /10A	74	67.00 IN	74.00 IN	83.00 IN
A-SCD-D-080P-SS-N-050	5 (3.7)	1.5 NPT	0.5 NPT	80	15.2	14.1	15.2 / 14.1	1844 LBS	34.4 (3.20)	12,722	36A / 31A/ 16A	74	67.00 IN	74.00 IN	83.00 IN
A-SCD-D-080P-SS-N-075	7.5 (5.6)	1.5 NPT	0.5 NPT	80	24	25.2	24 / 25.2	1934 LBS	34.4 (3.20)	19,083	49A/ 43A / 22A	74	67.00 IN	74.00 IN	83.00 IN