




America

CERTIFICATE

No. U8 098792 0003 Rev. 00

Holder of Certificate: **Amico Accessories Inc.**
 122 East Beaver Creek
 Richmond Hill ON L4B 1G6
 CANADA

Certification Mark:



Product: **General Medical Devices**
Portable Computer Cart / Mobile Workstation

This product was voluntarily tested to the relevant safety requirements referenced on this certificate. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL and a Standards Council of Canada accredited Certification body.

Test report no.: 7169009106-000

Date, 2022-02-08

(Peter Keith)



America

CERTIFICATE

No. U8 098792 0003 Rev. 00

Model(s): HMF-A-BBBCDD-EF(GG)-(H)(I)JKL,
HMDABBBBCDEFGHHIJKL,
HLX-BBC-DE-FF,
HMC-ABCC-DDEE-FF,
HML-A-BBB-CD(EE)-(F)(G)HI

(See below for model nomenclature)

Brand Name: **Hummingbird**

Tested according to: ANSI/AAMI ES60601-1:2005/A1:2012-08
CAN/CSA-C22.2 No. 60601-1:2014

Parameters:

Rated Input Voltage: 100-240 Vac
Rated Input Frequency: 50/60 Hz
Rated Input Current: 6.5 A
Protection Class: Class I

Excludes Usability (clauses 12.2 and 15.1)



America

CERTIFICATE

No. U8 098792 0003 Rev. 00

ZERTIFIKAT ◆ CERTIFICATE ◆ 認證書 ◆ CERTIFICADO ◆ CERTIFIKAT ◆ CERTIFICAT

HMF-	A-	BBB	C	DD-	E	F	(GG)-	(H)	(I)	J	K	L
<p>Product line (Hummingbird Mobile Workstation - For LCD)</p> <p>A = configuration, one of: P, Powered cart; or N, Non-powered cart If the Cart is non powered (N) then keep (GG), (H) & (I) blank & No ELE/ELL option</p> <p>BBB = Lifting, one of: 70,70N Gas spring,0-9.9lbs(0-4.49kg); 120,120N Gas spring,10-19.9lbs(4.54-9.03kg); 165,165N Gas spring,20-29.99lbs(9.07-13.6kg) 210,210N Gasspring,30-39.99lbs(13.61-18.14kg); 255, 255N Gas spring, 40 - 49.99lbs(18.15 - 22.68kg); or ELE, Actuator (Electric lift) 50 - 70lbs(22.68 - 31.75kg). (only with powered cart) ELL, Actuator (Electric lift) 50 - 70lbs(22.68 - 31.75kg) - with longer column. (only with powered cart)</p> <p>C =LCD mount L, Single LCD mount or 0, no LCD mount</p> <p>DD =Spring Strength, one of: 25, 5.5-8lbs(2.5-3.6kg)monitors & peripherals; 35, 7.5-10lbs(3.4-4.5kg)monitors & peripherals; 50, 13-15.5lbs(5.8-7.0kg)monitors & peripherals; 60, 15.5-18lbs(7-8.1kg)monitors & peripherals; 75, 18 - 22lbs(8.1 - 10kg) or lower monitors & peripherals; FH, Over 22Lbs(10kg) monitors & peripherals no height adjustments; 00, No FIT post or spring; or XX, Custom (monitors & peripherals >22lbs (10kg) or fixed height.</p> <p>E =Cart Power type, One of; U, North America plug type, 110V 60Hz (inlet) G, British plug, 230V 50Hz (inlet),I, Australia plug, 230V 50Hz (inlet) F, Schuko plug (inlet),K, Danish plug (inlet), X, Custom plug (inlet) 0, No Inlet</p> <p>F =Cart power outlet plug type U, North America plug type (outlet)/W, Universal (outlet)/X, Custom Power Outlet/0, No Power Outlet</p> <p>(GG)= Battery type, One of; 19, Powered with one LiFeoPO4 battery; 38, Powered with two LiFePO4 batteries; 57, Powered with three LiFePO4 batteries; 42, Powered with one SLA battery; or 20, Powered with One LiFePO4 Battery X-HMW-BAT-PSB122; 40, Powered with two LiFePO4 Battery X-HMW-BAT-PSB122;23, Powered with One LiFePO4 Battery X-HMW-BAT-PSL122; 46, Powered with two LiFePO4 Battery X-HMW-BAT-PSL122; or 0 0, None, If the cart is non powered, keep this section as blank</p> <p>(H) = Control Panel, one of; M, Standard interface; L, Advanced interface; or If the cart is non powered, keep this section as blank</p> <p>(I) = Cold Boot, One of; R, CPU Reset Button; or 0, None If the cart is non powered, keep this section as blank</p> <p>J = Handle type B, Low profile Handle</p> <p>K = Keyboard Tray, one of; S, Side out Keyboard tray; E, Ergonomic keyboard tray; X, Custom Keyboard tray; 0, None</p> <p>L = Bucket type, one of; X, Custom Bucket; 0, None</p>												



America

CERTIFICATE

No. U8 098792 0003 Rev. 00

ZERTIFIKAT ◆ CERTIFICATE ◆ 認證書 ◆ CERTIFICADO ◆ CERTIFIKAT ◆ CERTIFICAT

HMD	A	BBB	C	D	E	F	G	HH	(I)	J	K	L
Product line (Hummingbird CTM – IT workstation)												
<p>A = configuration, one of: P, Powered cart; or</p> <p>BBB = Lifting, one of: 70, 70N Gas spring, 0-9.9lbs(0 4.49kg); 120, 120N Gas spring, 10-19.9lbs(4.54-9.03kg); 165, 165N Gas spring, 20-29.99lbs(9.07-13.6kg); 210, 210N Gasspring, 30-39.99lbs(13.61-18.14kg); 255, 255N or lower xxN Gas spring, 40 - 49.99lbs(18.15 - 22.68kg); ELE, Actuator (Electric lift) 50 - 70lbs(22.68 - 31.75kg); or ELL, Actuator (Electric lift) 50 - 70lbs(22.68 - 31.75kg) with longer column</p> <p>C = CTM column length, one of: A, CTM Column of 18" length.</p> <p>D = Display orientation, one of: 5, Left side Arm</p> <p>E = Display orientation, one of: 5, Right side Arm</p> <p>F = Cart Power type, One of: G, British plug, 230V 50Hz (inlet) I, Australia plug, 230V 50Hz (inlet), F, Schuko plug (inlet), K, Danish plug (inlet), X, Custom plug (inlet), U, North America plug type, 110V 60Hz (inlet), 0, No Inlet</p> <p>G = Cart power outlet plug type U, North America plug type (outlet), W, Universal (outlet), X, Custom Power Outlet, 0, No Power Outlet</p> <p>HH = Battery type, One of: 19, Powered with one LiFeoPO4 battery; 38, Powered with two LiFePO4 batteries; 57, Powered with three LiFePO4 batteries, 42, Powered with one SLA battery; 20, Powered with One LiFEPO4 Battery X-HMW-BAT-PSB122; 40, Powered with two LiFEPO4 Battery X-HMW-BAT-PSB122; 23, Powered with One LiFEPO4 Battery X-HMW-BAT-PSL122; 46, Powered with two LiFEPO4 Battery X-HMW-BAT-PSL122; or 00, None</p> <p>I = Control Panel, One of: M, Standard interface; L, Advanced interface;</p> <p>J = Cold Boot, One of: R, CPU Reset Button;</p> <p>K = Handle type B, Low profile Handle</p> <p>L = Keyboard Tray, one of; S, Side out Keyboard tray; E, Ergonomic keyboard tray; X, Custom Keyboard tray; or 0, None</p>												



America

CERTIFICATE

No. U8 098792 0003 Rev. 00

ZERTIFIKAT ◆ CERTIFICATE ◆ 認證書 ◆ CERTIFIKAT ◆ CERTIFICADO ◆ CERTIFICAT

HL	X-	BB	C-	D	E-	FF
Product line (Hummingbird Lite – IT workstation) X = Configuration, one of: 1, Laptop; or 2, Laptop with Security Bracket, or 3, LCD/Display BB = Pneumatic Lift, one of: 16, 165N Gas spring, 3.5 lbs - 13.4 lbs (1.59 kg - 6.1 kg); 21, 210N Gas spring, 17 lbs - 22.5 lbs (10.21 kg - 17.71 kg); 25, 255N Gas spring, 13.5 lbs – 27 lbs (6.12 kg - 12.25 kg) C = LCD/Display Strength, one of: 3, 35N Gas spring, 7.5-10 lbs. (3.4-4.5 kg); 5, 50N Gas spring, 13-15.5 lbs. (5.8-7 kg); 6, 60N Gas spring, 15.5-18 lbs. (5.8-7 kg); 7, 75N Gas spring, 18-22 lbs. (8.1-19 kg); T, Tablet 5.5-8 lbs. (2.5-3.6 kg) 25N Gas spring (35N + 2.5 lbs Counter Weight) L, LCD/Display with Small form Factors, 5.5-15.5 lbs. (2.5-7kg) CPU 25N with spare 50N gas spring A, All In One 15.5-22 lbs. (5.8-19 kg) 60N with spare 75N gas spring FH, Fixed Height (blank), Laptop option (X = 1 or 2) D = Caster type, one of: D, 1 Directional Lock; M, 2 Total Locking; or T, 2 Total Locking (Tente). E = Mouse/Keyboard, one of: M, Standard Mouse Tray; S, Sliding Keyboard Tray with Mouse Tray); (For X = 3 ONLY) 0, None FF = RFID combinations, one of: XX, RFID combinations 00, None (blank), Laptop option (X = 1 or 2)						



America

CERTIFICATE

No. U8 098792 0003 Rev. 00

ZERTIFIKAT ◆ CERTIFICATE ◆ 認證證書 ◆ СЕРТИФИКАТ ◆ CERTIFICADO ◆ CERTIFICAT

HMC Product line ("CTM Cart – IT Work station")	A	B	CC-	DD	EE-	FF
<p>A =Base type, one of: L, Non-Powered with sheet metal base; or N, Non-Powered with casted base.</p> <p>B =Caster type, one of: 1, Two Tente total lock casters; 2, Two Dersheng total lock casters 3, One directional & one total lock casters (Tente); 4, One direction & one total lock caster (Dersheng); 5, Four Tente total lock casters; or 6, Four Dersheng total lock casters.</p> <p>CC = Extension length & color, one of: 3W, 48in overall height 5W, 60in overall height</p> <p>DD =Display Mount type, One of: L4, Flush LCD mount; LF, Ultra Flush LCD mount; LP, Laptop work surface; P5, PM5-MON; 40, SSM45 with 5" shelf & with 30deg tilt limit ; 42, SSM45 with 6" drop on Shelf; 43, SSM45 with Philips shelf; 44, SSM45 with accessory head; 47, SSM45 with 5" shelf & with slide in VESA plate; 48, SSM45 with 3" shelf & with 30deg tilt limit; 49, SSM45 with 3" shelf; 4W, SSM45 with welch Allyn Head; 70, 6" extension with 5" shelf w/30deg tilt limit; 73, 6" extension with Philip shelf; or 00, none.</p> <p>EE = Keyboard Mount type, One of: KF, Ultra flush keyboard; MF, Mini Keyboard; or 00, None.</p> <p>FF = Power options plug type & no. of plugs, one of: B4, North America with 4 plugs 15" cord length; G4, British plug I4, Australia plug F4, Schuko plug K4, Danish plug XX, Custom (other countries); or 00, None.</p>						



America

CERTIFICATE

No. U8 098792 0003 Rev. 00

ZERTIFIKAT ◆ CERTIFICATE ◆ CERTIFICADO ◆ CERTIFIKAT ◆ 認證書 ◆ CERTIFICATE ◆ ZERTIFIKAT

HML- Product line ("Hummingbird Mobile Workstation - For Laptop")	A-	BBB-	C	D	(EE)-	(F)	G	H	I
<p>A =Configuration, one of:P, Powered cart; or N, Non-powered cart If the Cart is non powered (N) then keep (GG), (H) & (I) blank & No ELE/ELL option</p> <p>BBB = Lifting, one of: 70,70N Gas spring,0- 9.9lbs(0-4.49kg); 120,120N Gas spring,10-19.9lbs(4.54-9.03kg); 165,165N Gas spring,20-29.99lbs(9.07-13.6kg) 210,210N Gasspring,30-39.99lbs(13.61- 18.14kg); 255, 255N Gas spring, 40 - 49.99lbs(18.15 - 22.68kg); or ELE, Actuator (Electric lift) 50 - 70lbs(22.68 - 31.75kg). (only with powered cart) ELL, Actuator (Electric lift) 50 - 70lbs(22.68 - 31.75kg) - with longer column. (only with powered cart)</p> <p>C = Cart Power type, One of; U, North America plug type, 110V 60Hz (inlet) G, British plug, 230V 50Hz (inlet) I, Australia plug, 230V 50Hz (inlet) F, Schuko plug (inlet) K, Danish plug (inlet) X, Custom plug (inlet) 0, No Inlet</p> <p>D = Cart power outlet plug type: U, North America plug type (outlet) W, Universal (outlet) X, Custom Power Outlet 0, No Power Outlet</p> <p>(EE) = Battery type, One of; 19, Powered with one LiFeoPO4 battery, 38, Powered with two LiFePO4 batteries; 57, Powered with three LiFePO4 batteries; 42, Powered with one SLA battery; 20, Powered with One LiFePO4 Battery X- HMW-BAT PSB122; 40, Powered with two LiFePO4 Battery X-HMW-BAT-PSB122; 23, Powered with One LiFePO4 Battery X-HMW-BAT-PSL122; 46, Powered with two LiFePO4 Battery X-HMW-BAT-PSL122; or 00, None , If the cart is non powered, keep this section as blank</p> <p>(F) = Control Panel, One of; M, Standard interface; L, Advanced interface; or If the cart is non powered, keep this section as blank</p> <p>(G)= Cold Boot, One of; R, CPU Reset Button; or 0, None If the cart is non powered, keep this section as blank</p> <p>H = Handle type B, Low profile Handle</p> <p>I =Keyboard Tray, one of; S, Side out Keyboard tray; E, Ergonomic keyboard tray; X, Custom Keyboard tray OR 0, None</p>									