

NULOOK SERENITY SERIES
SURFACE MOUNTED PANEL WALL
FULL LENGTH GENERAL CARE
(M/N: SFWXXXX-CPX)

FOR PRESENTATION
PURPOSES ONLY
Sep-17-2025
NOT FOR
CONSTRUCTION USE



TYPE:	QUANTITY:
PATIENT ROOM	00


IF MIRRORED UNITS ARE REQUIRED,
THE CUSTOMER MUST INDICATE THE
QUANTITY BELOW
QTY MIRRORED: _____

SYSTEM DETAILS

SYM	QTY	DESCRIPTION
PS	1	PANEL WALL
O	2	AMICO GAS OUTLET, OXYGEN
A	2	AMICO GAS OUTLET, MED AIR
V	2	AMICO GAS OUTLET, VACUUM
E2	2	RECEPTACLE - DUPLEX RED
BR	1	RECEPTACLE - BED IVORY
DA	2	PROVISION - DATA 1G
NC	1	PROVISION - NURSE CALL 3G
SLM	2	VACUUM SLIDE

IMPORTANT NOTES:

- STANDARD AVAILABLE WIDTHS OF SIDE PANELS: 19", 23.5", 32", 36", 42", 47.5" and 59".
- MAINTAIN A MINIMUM SPACING OF 3" C/C BETWEEN LOW VOLTAGE DEVICES, POWER OUTLETS, AND SWITCHES.
- MAINTAIN A MINIMUM SPACING OF 5" C/C BETWEEN GAS OUTLETS.
- NURSE CALL DEVICES MUST MAINTAIN A MINIMUM SPACING OF 6" C/C FROM OTHER SERVICES DUE TO THE SIZE OF THE DEVICE FACE



FOR REFERENCE ONLY

IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS
CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE

DATE

PHONE NO.

HOSPITAL	NURSE CALL	MFGR. : TBC	MODEL No. : TBC	 www.amico.com	85 Fulton Way Richmond Hill, Ontario L4B 2N4, CANADA Toll-Free: 1-877-462-6436(T) 1-866-440-4986(F) Tel: (905) 764-0800 Fax: (905) 764-0862	QUOTE No. REV. BY: RJUSAYAN PROJECT PM: DATE: Sep-17-2025 REV. NO:
LOCATION	GAS OUTLETS	MFGR. : AMICO	TYPE CONNECTION : TBC			
UNIT AS SHOWN - N.T.S.	FINISH	TBC	CEILING HEIGHT			

NULOOK SERENITY SERIES
SURFACE MOUNTED PANEL WALL
FULL LENGTH GENERAL CARE
(M/N: SFWXXXX-CPX)

FOR PRESENTATION
PURPOSES ONLY
Sep-17-2025
NOT FOR
CONSTRUCTION USE



SINGLE SIDED
UNIT DETAIL



SINGLE SIDED AS SHOWN SIDE

IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS
CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE

DATE

PHONE NO.

HOSPITAL	HOSPITAL	NURSE CALL	MFGR. : TBC	MODEL No. : TBC	 www.amico.com	85 Fulton Way Richmond Hill, Ontario L4B 2N4, CANADA Toll-Free: 1-877-462-6436(T) 1-866-440-4986(F) Tel: (905) 764-0800 Fax: (905) 764-0862	QUOTE No. REV. BY: RJUSAYAN PROJECT PM: DATE: Sep-17-2025 REV. NO:
LOCATION	LOCATION	GAS OUTLETS	MFGR. : AMICO	TYPE CONNECTION : TBC			
UNIT AS SHOWN - N.T.S.		FINISH	TBC	CEILING HEIGHT	TBC		