

REGAL SERIES
SURFACE MOUNTED FLATWALL
FULL LENGTH GENERAL CARE
(M/N: SFW-XXSR)



Install
Instructions

FOR PRESENTATION
PURPOSES ONLY
Sep-17-2025
NOT FOR
CONSTRUCTION USE

DRAWING # 00

TYPE:	QUANTITY:
PATIENT ROOM	00

IF MIRRORED UNITS ARE REQUIRED, THE
CUSTOMER MUST INDICATE THE QUANTITY
BELOW


QTY MIRRORED: _____

SYSTEM DETAILS

SYM	QTY	DESCRIPTION
FW	1	FLATWALL, XX" WIDTH
O	1	AMICO GAS OUTLET, OXYGEN
A	1	AMICO GAS OUTLET, MED AIR
V	1	AMICO GAS OUTLET, VACUUM
E2	2	RECEPTACLE - DUPLEX RED
NC	1	PROVISION - NURSE CALL 3G
DA	1	PROVISION - DATA 1G
SLM	1	VACUUM SLIDE

IMPORTANT NOTES:

- STANDARD AVAILABLE WIDTHS: 19", 23.5", 32", 36", 42", 47.5" and 59".
- MAINTAIN A MINIMUM SPACING OF 3" C/C BETWEEN LOW VOLTAGE DEVICES, POWER OUTLETS, AND SWITCHES.
- MAINTAIN A MINIMUM SPACING OF 5" C/C BETWEEN GAS OUTLETS.
- NURSE CALL DEVICES MUST MAINTAIN A MINIMUM SPACING OF 6" C/C FROM OTHER SERVICES DUE TO THE SIZE OF THE DEVICE FACE



FOR REFERENCE ONLY

IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS
CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE

DATE

PHONE NO.


HOSPITAL	HOSPITAL	NURSE CALL	MFGR. : TBC	MODEL No. : TBC	 www.amico.com	85 Fulton Way Richmond Hill, Ontario L4B 2N4, CANADA Toll-Free: 1-877-462-6436(T) 1-866-440-4986(F) Tel: (905) 764-0800 Fax: (905) 764-0862	QUOTE No. REV. BY: RJUSAYAN PROJECT PM: DATE: Sep-17-2025 REV. NO:
LOCATION	LOCATION	GAS OUTLETS	MFGR. : AMICO	TYPE CONNECTION : TBC			
UNIT AS SHOWN - N.T.S.		FINISH	TBC	CEILING HEIGHT			

REGAL SERIES
SURFACE MOUNTED FLATWALL
FULL LENGTH GENERAL CARE
(M/N: SFW-XXSR)

FOR PRESENTATION
PURPOSES ONLY
Sep-17-2025
NOT FOR
CONSTRUCTION USE



SINGLE SIDED
UNIT DETAIL



SINGLE SIDED AS SHOWN SIDE

IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS
CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE

DATE

PHONE NO.

HOSPITAL	HOSPITAL	NURSE CALL	MFGR. : TBC	MODEL No. : TBC	 www.amico.com	85 Fulton Way Richmond Hill, Ontario L4B 2N4, CANADA Toll-Free: 1-877-462-6436(T) 1-866-440-4986(F) Tel: (905) 764-0800 Fax: (905) 764-0862	QUOTE No. REV. BY: RJUSAYAN PROJECT PM: DATE: Sep-17-2025 REV. NO:
LOCATION	LOCATION	GAS OUTLETS	MFGR. : AMICO	TYPE CONNECTION : TBC			
UNIT AS SHOWN - N.T.S.		FINISH	TBC	CEILING HEIGHT			
				TBC			