




Features

- 60" reach for optimal patient experience
- Ergonomic design for reduced clinician fatigue
- VESA 75/100 mounting pattern supports the wide range of devices sizing up to 24"
- Built-in cable management for hassle free operation
- 2 extended arms provided to optimal reach and articulation
- Maximum Load: 2.21 to 14.33 lbs (1 to 6.5 kg)
- Designed to meet infection control and prevention standards using non-porous, easy-to-clean materials that are safe for use with hospital grade cleaning solutions

Specifications

Cable Routing

- USB-A, USB-C, PoE++ applications, Coaxial, Cat6, any cable with touch-safe voltage supply



WARNING: DO NOT ROUTE AC POWER CORD THROUGH THE ARMS.

Display Articulation

- Hinge Swivel: +/- 185°
- Device Mount Rotation: +/- 5°
- Hinge Tilt: -60° to +20°

Display Support

- Mounting Pattern: 100 mm x 100 mm / 75 mm x 75 mm
- Mount: Single Device
- Display Type: Single Monitor
- Panel Size Support: 15" to 24"
- Support Weight Range: 2.21 to 14.33 lbs (1 to 6.5 kg)
- Arm 1 Swivel: +/- 90°

Height Adjustment

- Height Adjustment Range: 57.76" (146.7 cm)
- Tension Torque Adjust: Yes
- Arm 1 Tilt: 90°
- Arm 2 Tilt: 82.5°

Key Dimensions


Outside Dimensions (in extended position):

- Width: 7.87" (20.0 cm)
- Depth: 73.74" (187.3 cm)
- Height: 82.91" (210.6 cm)


Outside Dimensions (in folding position):

- Width: 7.87" (20.0 cm)
- Depth: 13.62" (34.6 cm)
- Height: 54.1" (137.4 cm)
- Physical Weight: 28.44 lbs (12.9 kg)

Figure 1: Long Reach Arm Articulation and Dimensions



Please visit www.amico.com for all technical specifications, installation instructions, detailed AutoCAD drawings, and more Mounting Solutions.



Distributed By: