

Hummingbird Accessories Catalog



Hummingbird Cart



Perch Accessories Rail 3

- Recommended for the powered Hummingbird cart with a monitor weighing less than 21 lbs (9.5 kg) or non-powered Hummingbird cart with a monitor weighing a maximum of 9 lbs (4.1 kg).
- The maximum load is 5 lbs (2.3 kg) with one Perch Rail accessory on each side or 10 lbs (4.5 kg) of one accessory in the middle only.

VESA Mount 4

- Compatible with Tablet/LCD/AIO carts only. Recommended to use one VESA accessory per Hummingbird cart.
- The maximum load depends on the weight of the secondary lifting capacity minus the weight of the Tablet/LCD/AIO monitor. If the LCD mount is at the fixed height then the maximum load depends on the cart lifting capacity.

Worktop & Main Bin 5

- Recommended for the Tablet/LCD/AIO carts only.
- The maximum load depends on the weight of the secondary cart lifting capacity minus the weight of the Tablet/LCD/AIO monitor or 30 lbs (13.6 kg), whichever is lower.

Bump Accessories Rail 6

- The maximum load is 5 lbs (2.3 kg) with one Bump Rail accessory on each side or one 10 lbs (4.5 kg) accessory in the middle only.
- May increase footprint depending on the accessory type.

Column 7

- The maximum load depends on the lifting capacity of the cart if there is an Ergonomic Rail added. Compatibility also depends on the height of the accessories already selected.

Index 10

- Hummingbird Carts & Compatibility

Compatible Cart Type

Mounting Platforms	Hummingbird Lite HL1, HL2, & HL3	CTM Cart HMC	Hummingbird Worktop HMD	Hummingbird LCD HMF	Hummingbird Laptop HMLP
Perch Rail				●	
VESA Mount	○	○○○		●	
Worktop & Main Bin	○○		●	●	●
Bump Rail		●	●	●	●
Column	●	●	●	●	●

● Available

○ Available to HL3 only

○○ Available with CTM Cart for LCD/Tablet/AIO only

○○○ Available to with Worktop accessories only

Perch Accessory Rail

Compatible with Hummingbird LCD/Tablet/AIO Cart
















HMA-PR-A-18

18" (45.7 cm) Perch Accessory Rail for Hummingbird LCD/Tablet/AIO cart on top of the monitor post. Populate with one or two Perch accessories for easy access. Recommended to work with the Powered Hummingbird cart with a monitor that has a maximum weight of 21 lbs (9.5 kg), or a non-powered Hummingbird cart with a monitor that has a maximum weight of 9 lbs (4 kg).

For the Perch Accessory Rail & Accessories by Adapter Orientation

Compatible with Hummingbird LCD Cart

Hummingbird Perch accessories mount to the Perch Accessory Rail (HMA-PR-A-18). Can be configured in various orientations to fit your unique workflow. Each Perch Rail is available with up to 2 accessories. The maximum load for the 5" x 5" x 4" (12.7 cm x 12.7 cm x 10 cm) basket is 4.5 lbs (2 kg), rest at 1.1 lbs (0.5 kg).

Adapter Orientation		Accessories	
	HMA-(X)-(FTA) On the Perch.		HMA-050504-FTA For small supplies that fit into 5" x 5" x 4" (12.7 cm x 12.7 cm x 10.1 cm) dimensions with a maximum load of 4.5 lbs (2 kg).
			HMA-BARCD-FTA Ideal for holding a tethered barcode scanning gun. Barcode scanner is for illustration purposes only, not included in price.
			HMA-IVHOOK2-FTA Designed to hold IV bags.
	HMA-(X)-(LA9) Perch with 9" drop, 90° turn capabilities, makes it easy to access from the side.		HMA-050504-LA9 For small supplies that fit into 5" x 5" x 4" (12.7 cm x 12.7 cm x 10.1 cm) dimensions with a maximum load of 4.5 lbs (2 kg).
			HMA-BARCD-LA9 Ideal for holding a tethered barcode scanning gun. Barcode scanner is for illustration purposes only, not included in price.
			HMA-BARCDHLD2-LA9 Designed to fit a wireless tethered barcode scanning gun.
	HMA-(X)-(FD9) Perch with 9" drop.		HMA-050504-FD9 For small supplies that fit into 5" x 5" x 4" (12.7 cm x 12.7 cm x 10.1 cm) dimensions with a maximum load of 4.5 lbs (2 kg).
			HMA-BARCD-FD9 Ideal for holding a tethered barcode scanning gun. Barcode scanner is for illustration purposes only, not included in price.
	HMA-(X)-(LSA) Perch with extended reach capabilities.		HMA-BARCD-LSA Ideal for holding a tethered barcode scanning gun right on the side of the monitor while facing the user. Barcode scanner is for illustration purposes only, not included in price.

Accessories for the VESA Mount

Compatible with Hummingbird LCD/Lite (HL3)/CTM Cart



MON-LCD-CPUMNT

Mounts a CPU with dimensions of 10" x 10" x 3" (25.4 cm x 25.4 cm x 7.62 cm) directly behind an LCD/flat panel display.

LCD monitor and CPU are for illustration purposes only, not included in price.



MON-LCD-CPUMNT-SM

Mounts a VESA compatible CPU with dimensions of 10" x 10" x 1.75" (25.4 cm x 25.4 cm x 4.45 cm) directly behind an LCD/flat panel display.

LCD monitor and CPU are for illustration purposes only, not included in price.



HMA-BARCDHLD2-LCD

Mounts a wireless barcode scanner dock from behind to the side of an LCD/flat panel display.

LCD monitor & barcode scanner are for illustration purposes only, not included in price.



MON-LCD-BARCD

Mounts a tethered barcode scanner gun from behind to the side of an LCD/flat panel display.

LCD monitor & barcode scanner are for illustration purposes only, not included in price.

























Accessories for the Worktop & Main Bin

Compatible with Hummingbird LCD/Laptop/Worktop Carts

	HMA-WT-LEDMNT Hummingbird interface extension 12 Volt LED light for the worktop of the Hummingbird cart.		HMA-WT-USB-4HUB 4-hub USB splitter for the Hummingbird cart.
	HMA-USBEXT-3 Extended USB cord.		HMA-WRISTREST-CVR Replaceable wrist rest cover for the Hummingbird cart.
	MONWRISTPAD Antimicrobial comfort wrist rest.		MONWRISTPAD-S Small antimicrobial comfort wrist rest.
	HMA-MSBIN-KIT-E Mouse holder for the Hummingbird cart. Mouse is for illustration purposes, not included in price.		HMA-SURFLABEL-(XXXX) Personalized surface label for the Hummingbird worktop for enhancing the branding of your facility. Generic surface label is for illustration purposes only.
	HMA-RFID-(XXXX) RFID system for the Hummingbird cart. Generic SKU for quoting only. Please call for details.		


Accessories for the Worktop & Main Bin

Compatible with Hummingbird Lite Cart

	HMA-UNI-BARCD-PWRH Tethered barcode scanner gun holder, placed at the bottom of the worktop for the Hummingbird Lite for Laptop (integrated secure bracket model only). Barcode scanner is for illustration purposes only, not included in price.		MON-UNI-BARCD Tethered barcode scanner gun holder, placed at the exiting power cord holder is designed for the Hummingbird Lite LCD cart (HL3), or the Laptop cart with no secure bracket (HL1).																		
	X-HLT-PWR-HOOK The power cord holder keeps the cables organized. Designed for the Hummingbird Lite laptop cart (HL1) with no integrated secure bracket model, or the Hummingbird Lite LCD cart (HL3).	 <small>* Generic part number for quoting purposes, part number may vary.</small>	HMA-WPWRBAR-(X) International power bar with C14 coiled power cord, inlet plug to match local plug types (perfect for the Hummingbird Lite cart). *																		
	X-HLT-WTP-ASY-___ RFID reader for the Humminbird Lite LCD/ Tablet/AIO cart to be placed at the base of the IT post in the worktop.	<table><tr><th>Plug</th><th>Plug Type</th><th>Part Number</th></tr><tr><td></td><td>Type B (US)</td><td>HMA-WPWRBAR-B</td></tr><tr><td></td><td>Type G (UK)</td><td>HMA-WPWRBAR-G</td></tr><tr><td></td><td>Type I (Australia)</td><td>HMA-WPWRBAR-I</td></tr><tr><td></td><td>Type E, F (Schuko)</td><td>HMA-WPWRBAR-F</td></tr><tr><td></td><td>Type K (Denmark)</td><td>HMA-WPWRBAR-K</td></tr></table>		Plug	Plug Type	Part Number		Type B (US)	HMA-WPWRBAR-B		Type G (UK)	HMA-WPWRBAR-G		Type I (Australia)	HMA-WPWRBAR-I		Type E, F (Schuko)	HMA-WPWRBAR-F		Type K (Denmark)	HMA-WPWRBAR-K
Plug	Plug Type	Part Number																			
	Type B (US)	HMA-WPWRBAR-B																			
	Type G (UK)	HMA-WPWRBAR-G																			
	Type I (Australia)	HMA-WPWRBAR-I																			
	Type E, F (Schuko)	HMA-WPWRBAR-F																			
	Type K (Denmark)	HMA-WPWRBAR-K																			
	MONWRISTPAD Antimicrobial comfort wrist rest.																				
	MONWRISTPAD-S Small antimicrobial comfort wrist rest.																				

Accessories for the Bump

Compatible with Hummingbird LCD/Laptop/Worktop Cart



HMA-100609-BMP

10" x 6" x 9" (25.4 cm x 15.2 cm x 22.8 cm) basket on the Bump for the Hummingbird cart. (Ref: old part # HMA-BMP-1-0-E-0)



HMA-CHARTHLD-BMP

Standard chart holder for file folders up to 12" (30.4 cm) wide on the Bump for the Hummingbird cart. (Ref: old part # PROJ-0724-1)
Folder is for illustration purpose, not included in price.

Bump Accessory Rail

Compatible with Hummingbird LCD/Laptop/Worktop Cart



HMA-A18-BMP

Standard 18" (45.7 cm) Bump Accessory Rail on the Hummingbird cart. Maximum weight is 5 lbs (2.24 kg) static load on each side of the rail. Heavier accessories can be mounted at the center of the rail if it is the only accessory on the rail (up to 10 lbs [4.5 kg] static load).

Accessories for the Bump Accessory Rail

Compatible with Hummingbird LCD/Laptop/Worktop Cart



100609-F-BB *

Can fit in 10" x 6" x 9" (25.4 cm x 15.2 cm x 22.8 cm) dimensions, with a maximum load of 2 lbs (0.90 kg).



190503-F-BT *

Can fit in 19" x 5" x 3" (48.4 cm x 12.7 cm x 7.6 cm) dimensions with a maximum load of 2 lbs (0.90 kg).



050504-F-BB

Can fit in 5" x 5" x 4" (12.7 cm x 12.7 cm x 10.1 cm) dimensions with a maximum load of 1 lb (0.4 kg).



AM-CUFBSK-F

For small objects like a cuff and inflation bulb with a maximum load of 1 lb (0.4 kg).

Blood pressure cuff is for illustration purpose, not included in price.



CHARTHLD-F

For file folders up to 12" (30.4 cm) wide with a maximum load of 1.5 lbs (0.6 kg).

Folder is for illustration purpose, not included in price.



GLBOXHLD-1-F *

Single glove box holder 10" x 5" x 4" (25.4 cm x 12.7 cm x 10.1 cm) with a maximum weight of 2 lbs (0.9 kg).

Glove box is for illustration purpose, not included in price.



WIPEHLD-N-ADJ-F

For sanitizing wipe bottles (Metrex, PDI SKU-111100, Q89072, and Q8504 disinfection towelettes - 65,160 pieces) with a maximum load of 2.5 lbs (1.1 kg).

Sanitizer wipes dispenser is for illustration purpose, not included in price.



SANDSP-UNI-F *

For touchless soap/hand sanitizer dispensers weighing a maximum of 10 lbs (4.5 kg) (eg. Kimberly Clark 92145, GOJO 192004, Purell H2288, and Ecolab 61439813).



PWRHLD-F *

Holds a power adapter with a maximum load of 5 lbs (2.25 kg).

Power adapter is for illustration purpose, not included in price.



BARCD-F

Holds a tethered barcode scanning gun less than, or equal to 1 lb (0.4 kg).

Barcode scanner is for illustration purpose, not included in price.



BARCDHLD2-F

Holds wireless barcode scanner gun docking station with a maximum load of 2 lbs (0.9 kg).








Wireless barcode scanner and dock are for illustration purpose, not included in price.

* Recommended to be the only Bump accessory on a non-powered Hummingbird LCD cart with an LCD monitor weighing >= 10 lbs (4.5 kg).

Accessories for the Column

Compatible with Hummingbird LCD/Laptop/Lite/Worktop Cart

Mounting to the Front of the Column (Facing the User)









	<p>HMA-CPUMNT1-C</p> <p>Designed to house CPUs with a depth of 1.5" (4 cm) to 3" (7.5 cm). With a maximum load of 20 lbs (9.1 kg), or depends on the cart lifting capacity, whichever is lower.</p>		<p>HMA-CPUMNTA-C</p> <p>Designed to house CPUs with a depth of 1.5" (4 cm) to 7" (18 cm). With a maximum load of 20 lbs (9.1 kg), or depends on the cart lifting capacity, whichever is lower.</p>
	<p>HMA-101004-C</p> <p>10" x 10" x 4" (25.4 cm x 25.4 cm x 10.1 cm) medium basket, great for miscellaneous storage and organizing supplies. With a maximum load of 10 lbs (4.5 kg), or depends on the cart lifting capacity, whichever is lower.</p>		<p>HMA-100609-C</p> <p>10" x 6" x 9" (25.4 cm x 15.2 cm x 22.8 cm) medium basket that works well with small items can also work as a convenient waste receptacle. With a maximum load of 10 lbs (4.5 kg), or depends on the cart lifting capacity, whichever is lower.</p>
	<p>HMA-CHARTHLD-C</p> <p>Efficient patient chart organizer fits file folders up to 12" (30.5 cm) wide. With a maximum load of 1.5 lbs (0.7 kg), or depends on the cart lifting capacity, whichever is lower.</p> <p>Folder is for illustration purpose, not included in price.</p>		<p>HMA-SHELF1220</p> <p>12" x 20" (30.5 cm x 50.8 cm) shelf. With a maximum load of 5 lbs (2.3 kg), or depends on the cart lifting capacity, whichever is lower.</p>
	<p>HMA-ERG-RAIL-ACSS</p> <p>Hummingbird 18" (45.7 cm) Ergonomic Accessory Rail on the Column. With a maximum weight of 5 lbs (2.3 kg) static load on each side of the rail. Heavier accessories, up to 10 lbs (4.5 kg) static load, if being mounted at the center and is the only accessory.</p>		



Accessories for the Ergonomic Rail

Compatible with Hummingbird LCD/Laptop/Lite/Worktop Cart

Mounting to the Ergonomic Accessories Rail (Facing the User)

 <p>100609-F-BB</p> <p>For small supplies that fit in 10" x 6" x 9" (25.4 cm x 15.2 cm x 22.8 cm) dimensions, with a maximum load of 2 lbs (0.9 kg).</p>	 <p>101004-F-BB</p> <p>For small supplies that fit in 10" x 10" x 4" (25.4 cm x 25.4 cm x 10.1 cm) dimensions, with a maximum load of 2 lbs (0.9 kg).</p>
 <p>190503-F-BB</p> <p>For small supplies that fit in 19" x 5" x 3" (48.4 cm x 12.7 cm x 7.6 cm) dimensions, with a maximum load of 2 lbs (0.9 kg).</p>	 <p>050504-F-BB</p> <p>For small supplies that fit in 5" x 5" x 4" (12.7 cm x 12.7 cm x 10.1 cm) dimensions, with a maximum load of 1 lb (0.4 kg).</p>
 <p>HMA-050504-FE</p> <p>For small supplies that fit in 5" x 5" x 4" (12.7 cm x 12.7 cm x 10.1 cm) dimensions, with a maximum load of 1 lb (0.4 kg). Placed under the main bin extended out to the side to make it easy to access.</p>	 <p>AM-CUFBSK-F</p> <p>For small objects like cuff and inflation bulb, with a maximum load of 1 lb (0.4 kg). Cuff is for illustration purposes, not included in price.</p>
 <p>HMA-AM-CUFBSK-FE</p> <p>For small objects like cuff and inflation bulb, with a maximum load of 1 lb (0.4 kg).</p>	 <p>CHARTHLD-F</p> <p>For file folders up to 12" (30.5 cm) wide, with a maximum load of 1.5 lbs (0.7 kg). Folder is for illustration purposes, not included in price.</p>
 <p>HMA-CHARTHLD-FE</p> <p>For file folders up to 12" (30.5 cm) wide, with a maximum load of 1.5 lbs (0.7 kg). Placed under the main bin extended out to the side to make it easy to access. Folder is for illustration purposes, not included in price.</p>	 <p>GLBOXHLD-1-F</p> <p>For a glove box within 10" x 5" x 4" (25.4 cm x 12.7 cm x 10.1 cm), with a maximum load of 2 lbs (0.9 kg). Glove box is for illustration purposes, not included in price.</p>
 <p>SHARPSHLD-F</p> <p>Universal sharps holder. Placed under the main bin facing the user. Sharps container is for illustration purposes, not included in price.</p>	 <p>DA-SHARPSHLD-F</p> <p>Designed to work with Daniels SharpSmart Sharps Holder (product code: 4005007000).</p>
 <p>WIPEHLD-N-ADJ-F</p> <p>For sanitizing wipe bottles (Metrex, PDI SKU-111100, Q89072, and Q8504 disinfection towelettes - 65,160 pieces) with a maximum load of 2.5 lbs (1.1 kg). Sanitizing wipe bottle is for illustration purpose, not included in price.</p>	 <p>HMA-BARCDHLD2-FE</p> <p>For wireless barcode scanner gun docking station, with a maximum load of 2 lbs (0.9 kg). Placed under the main bin extended out to the side to make it easy to access. Barcode scanner is for illustration purpose, not included in price.</p>
 <p>SANIHLD-F</p> <p>For mounting squeeze-top bottle hand sanitizer. Hand sanitizer is for illustration purpose, not included in price.</p>	 <p>HMA-ZB-QL220-FE</p> <p>Mounts a Zebra QL220 label printer. Placed under the main bin extended out to the side to make it easy to access.</p>

Accessories for the Column

Compatible with Hummingbird LCD/Laptop/Lite/Worktop Cart

Mounting to the Side of the Column (Left or Right)



HMA-050504-V-BB

Small basket at 5" x 5" x 4" (12.7 cm x 12.7 cm x 10.1 cm) vertical mounted for small supplies or acts as a trash container.



HMA-100609-V-BB

10" x 6" x 9" (25.4 cm x 15.2 cm x 22.8 cm) basket, vertical mounting for the Hummingbird cart Column, side mounting.



HMA-101004-V-BB

Shallow basket at 10" x 10" x 4" (25.4 cm x 25.4 cm x 10.1 cm) vertical mounted for small supplies or acts as a trash container.



CHARTHLD-V

Efficient patient chart organizer fits file folders up to 12" (30.5 cm) wide. Mount on to the Hummingbird cart Column side. Maximum load of 1.5 lbs (0.7 kg).

Folder is for illustration purpose, not included in price.



HMA-CUFBASK-V

Cuff basket for small objects like a cuff and inflation bulb.

Cuff and inflation bulb are for illustration purposes, not included in price.



HMA-SANIHL-D-V

Hand sanitizer holder (squeeze-top bottle) for the Hummingbird cart Column, side mounting.

Bottle of hand sanitizer is for illustration purpose, not included in price.



HMA-SHELF5-06FLD-V

6" x 6" (15.2 cm x 15.2 cm) small swivel folding shelf.



HMA-GLBOXHLD-1-VE

Glove box holder for the Hummingbird cart Column, side mounting.

Glove box is for illustration purposes, not included in price.



WIPEHLD-N-ADJ-V

Ideal for mounting a sanitizing wipe bottle (like Metrex & SKU 1111100, Q89072, and Q8504 Disinfection Towelletes - 65,160 pieces). Maximum load is 2.5 lbs (1.1 kg).

Sanitizing wipe bottle is for illustration purpose, not included in price.



DA-SHARPSHLD-V

Designed to work with Daniels Sharpsmart Sharps Holder (product code: 4005007000) and mounts to the Hummingbird Column, side mounting.



HMA-SHARPSHLD-V

Universal sharps holder to keep a small sharps container within reach.

Sharps container is for illustration purposes, not included in price.



HMA-ZB-QL220-V

Mounts a Zebra QL220 Printer to the Hummingbird cart Column side Rail.



BARCDHLD2-V

Ideal for the wireless barcode scanner gun docking station, with a maximum load of 2 lbs (0.9 kg).

Barcode scanner is for illustration purpose, not included in price.



HMA-WPWRBAR-(X)-VE

International power bar with C14 coiled power cord, inlet plug to match local plug types (perfect for the Hummingbird Lite cart). *

* Generic part number for quoting purposes, part number may vary.



PWRBAR4-V

Power bar with 4 outlets, hospital grade 120 VAC 60 Hz, and a 15 ft (4.6 m) cord.



HMA-PWRBAR4-VE

Power bar with 4 outlets, hospital grade 120 VAC 60Hz, and a 15 ft (4.6 m) cord. Side mounting with offset plate towards/away from the user.

Plug	Plug Type	Part Number
	Type B (US)	HMA-WPWRBAR-B-VE
	Type G (UK)	HMA-WPWRBAR-G-VE
	Type I (Australia)	HMA-WPWRBAR-I-VE
	Type E, F (Schuko)	HMA-WPWRBAR-F-VE
	Type K (Denmark)	HMA-WPWRBAR-K-VE

Hummingbird Carts & Compatibility

HL1

Hummingbird Lite Laptop Cart



HL2

Hummingbird Lite Laptop Cart
with Integrated Secure Bracket



HL3

Hummingbird Lite LCD/AIO/
Tablet Cart



HMC

Hummingbird CTM Cart



HMD

Hummingbird Worktop Cart



HMF

Hummingbird LCD/AIO/Tablet
Cart



HMF-P

Hummingbird Powered LCD/AIO/
Tablet Cart



HML

Hummingbird Laptop Cart



Accessories	Mounting Surface	HL1	HL2	HL3	HMC	HMD	HMF	HMF-P	HML	Part Number
RFID Reader	In Worktop			●						X-HLT-WTP-ASY-__
	In Main Bin					●	●	●	●	HMA-RFID-XXXX
Wrist Rest - Small	On Worktop	●	●	●		●	●	●	●	MONWRISTPAD-S
Wrist Rest	On Worktop	●	●	●		●	●	●	●	MONWRISTPAD
LED Light for Worktop	On Worktop	●	●	●		●	●	●	●	HMA-WT-LEDMNT
USB 4-hub Splitter	In Main Bin					●	●	●	●	HMA-WT-USB-4HUB
Extended USB Cords	In Main Bin					●	●	●	●	HMA-USBEXT-3
Replaceable Wrist Rest Cover	On Keyboarding Surface	●	●	●		●	●	●	●	HMA-WRISTREST-CVR
Mouse Holder	Under Keyboard Tray					●	●	●	●	HMA-MSBIN-KIT-E
Additional USB Cooling Fan Inside CPU	In Main Bin					●	●	●	●	HMA-COOLINGFAN
Personalized Worktop Surface Label	On Worktop					●	●	●	●	HMA-SURFLABEL-XXXX
18" Perch Accessory Rail	Top of FIT Post							●		HMA-PR-A-18
18" Bump Accessory Rail	Bump					●	●	●	●	HMA-A18-BMP
18" Ergonomic Rail	Column Direct Mount (front)	●	●	●		●	●	●	●	HMA-ERG-RAIL-ACSS
	CTM				●					LRS-018-M
Power Cord	Bottom of Worktop	●		●						X-HLT-PWR-HOOK
Tethered Barcode Scanning Gun	Bottom of Worktop + Power Cord Hook		●							HMA-UNI-BARCD-PWRH
	Power Cord Hook	●		●						MON-UNI-BARCD
	Perch Rail							●		HMA-BARCD-FTA
								●		HMA-BARCD-LA9
								●		HMA-BARCD-FD9
								●		HMA-BARCD-LSA
	Bump Rail					●	●	●	●	BARCD-F
	VESA Mount			●	●	●	●	●		MON-LCD-BARCD
Wireless Barcode Scanner Dock	Perch Rail							●		HMA-BARCDHLD2-LA9
	Column Side Mount (left/right)	●	●	●		●	●	●	●	BARCDHLD2-V
	Ergonomic Rail + Extension	●	●	●		●	●	●	●	HMA-BARCDHLD2-FE
	CTM				●					BARCDHLD2-M
	Bump Rail					●	●	●	●	BARCDHLD2-F
	VESA Mount			●		●	●	●		HMA-BARCDHLD2-LCD
Swivel Shelf (6" x 6")	Column Side Mount (left/right)	●	●	●		●	●	●	●	HMA-SHELFS-06FLD-V
Printer Shelf	Column Direct Mount (front)	●	●	●		●	●	●	●	HMA-SHELF1220
	CTM				●					MPN-0171-ST
IV Hook	Perch Rail							●		HMA-IVHOOK2-FTA
AHM Handle	CTM				●					HMC-HAND

Accessories	Mounting Surface	HL1	HL2	HL3	HMC	HMD	HMF	HMF-P	HML	Part Number
Chart Holder	Column Side Mount (left/right)	●	●	●		●	●	●	●	CHARTHLD-V
	Column Direct Mount (front)	●	●	●		●	●	●	●	HMA-CHARTHLD-C
	Ergonomic Rail	●	●	●		●	●	●	●	CHARTHLD-F
	Ergonomic Rail + Extension	●	●	●		●	●	●	●	HMA-CHARTHLD-FE
	CTM				●					CHARTHLD-MON
	Bump					●	●	●	●	HMA-CHARTHLD-BMP
	Bump Rail					●	●	●	●	CHARTHLD-F
Printer (e.g. Zebra QL220 Label Printer)	Column Side Mount (left/right)	●	●	●		●	●	●	●	HMA-ZB-QL220-V
	Ergonomic Rail + Extension	●	●	●	●	●	●	●	●	HMA-ZB-QL220-FE
	CTM				●					BARCDHLD2-MON + ZB-QL220- PLT-00
	VESA Mount			●		●	●	●		MPN-0036 -LW
Sharps Container	Column Side Mount (left/right)	●	●	●		●	●	●	●	HMA -SHARPSHLD-V
	Ergonomic Rail	●	●	●	●	●	●	●	●	SHARPSHLD-F
	CTM				●					SHARPSHLD-MON
Daniel Sharps Container	Column Side Mount (left/right)	●	●	●		●	●	●	●	DA-SHARPSHLD-V
	Ergonomic Rail	●	●	●	●	●	●	●	●	DA-SHARPSHLD-F
	CTM				●					DA-SHARPSHLD-MON
Glove Box Holder (single)	Column Side Mount + Offset (left/right)	●	●	●		●	●	●	●	HMA-GLBOXHLD-1-VE
	Ergonomic Rail	●	●	●	●	●	●	●	●	GLBOXHLD-1-F
	CTM				●					GLBOXHLD-1-MON
	Bump Rail					●	●	●	●	GLBOXHLD-1-F
Squeeze-top Hand Sanitizer Bottle	Column Side Mount (left/right)	●	●	●		●	●	●	●	HMA-SANIHLD-V
	Ergonomic Rail	●	●	●	●	●	●	●	●	SANIHLD-F
	CTM				●					SANIHLD-M
Touchless Hand Sanitizer/ Soap Dispenser	Bump Rail					●	●	●	●	SANDSP-UNI-F
Sanitizing Wipe Can Holder	Column Side Mount (left/right)	●	●	●		●	●	●	●	WIPEHLD-N-ADJ-V
	Ergonomic Rail	●	●	●	●	●	●	●	●	WIPEHLD-N-ADJ-F
	CTM				●					WIPEHLD-N-ADJ-MON
	Bump Rail					●	●	●	●	WIPEHLD-N-ADJ-F
Cuff Basket	Column Side Mount (left/right)	●	●	●		●	●	●	●	AM-CUFBSK-V
	Ergonomic Rail	●	●	●	●	●	●	●	●	AM-CUFBSK-F
	Ergonomic Rail + Extension	●	●	●	●	●	●	●	●	HMA-AM-CUFBSK-FE
	CTM				●					AM-CUFBSK-M
	Bump Rail					●	●	●	●	AM-CUFBSK-F



Accessories	Mounting Surface	HL1	HL2	HL3	HMC	HMD	HMF	HMF-P	HML	Part Number
5" x 5" x 4" Basket	Perch Rail							●		HMA-050504-FTA
	Perch Rail							●		HMA-050504-LA9
	Perch Rail							●		HMA-050504-FD9
	Column Side Mount (left/right)	●	●	●		●	●	●	●	050504-V-BB
	Ergonomic Rail	●	●	●	●	●	●	●	●	050504-F-BB
	Ergonomic Rail + Extension	●	●	●	●	●	●	●	●	HMA-050504-FE
	CTM				●					050504-MON-BB
	Bump Rail					●	●	●		050504-F-BB
10" x 6" x 9" Basket	Column Direct Mount (front)	●	●	●		●	●	●	●	HMA-100609-C
	Column Side Mount (left/right)	●	●	●		●	●	●	●	HMA-100609-V-BB
	Ergonomic Rail	●	●	●	●	●	●	●	●	100609-F-BB
	CTM				●					100609-MON-BB
	Bump					●	●		●	HMA-100609-BMP
	Bump Rail					●	●	●	●	100609-F-BB
10" x 10" x 4" Basket	Column Direct Mount (front)	●	●	●		●	●	●		HMA-101004-C
	Column Side Mount (left/right)	●	●	●		●	●	●		HMA-101004-V-BB
	Ergonomic Rail	●	●	●	●	●	●	●		101004-F-BB
	CTM				●					101004-MON-BB
19" x 5" x 3" Basket	Bump Rail					●	●	●	●	190503-F-BT
	Ergonomic Rail	●	●	●	●	●	●	●	●	190503-F-BT
International Power Bar	Column Side Mount + Offset (left/right)	●	●	●		●	●		●	HMA-WPWRBAR(X)-VE
	Bottom of Worktop	●	●	●						HMA-WPWRBAR(X)
Power Bar	Column Side Mount (left/right)	●	●	●		●	●	●	●	PWRBAR4-V
	Column Side Mount + Offset (left/right)	●	●	●		●	●	●	●	HMA-PWRBAR4-VE
	CTM				●					PWRBAR4-MON
CPU (1.5" - 3" depth)	Column Direct Mount (front)					●	●	●		HMA-CPUMNT1-C
	CTM				●					CPUMNT1-MON
CPU (1.5" - 7" depth)	Column Direct Mount (z)					●	●	●		HMA-CPUMNTA-C
	CTM				●					CPUMNTA-M
Slim CPU (10" x 10" x 3") + VESA Compatible Mounting Case/Sleeves	VESA Mount					●	●	●		MON-LCD-CPUMNT
Slim CPU (10" x 10" x 1.75") + VESA Compatible Mounting Case/Sleeves	VESA Mount					●	●	●		MON-LCD-CPUMNT-SM
Power Adapter Holder	Bump Rail					●	●	●		PWRHLD-F
	CTM				●					PWRHLD-MON
Power Adapter Holder (dual)	CTM				●					PWRHLD2-MON




**Telehealth
Camera Mount**



SmartDrawers



Registration Cart



Amico Accessories Inc. | 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada
Toll Free Tel: 1.877.264.2697 | Tel: 905.763.7778 | Fax: 905.763.8587
Email: info@amico-accessories.com | www.amico.com

Represented By: