

Part Number: RS5-420-PKG



Part Number: RS5-LCD

Features:

- Cord wrap, basket and handle included in package models
- Compatibility with Welch Allyn, GE and Philips mounts
- Variety of accessories available

Specifications:

- Up to 35 lbs (15.8 kg) load capacity based on model
- Overall height 42" (106.7 cm)
- 3" (7.6 cm) locking casters
- Sturdy 24" base (61 cm)



Counterweights are available.

Monitor Compatibility

PATIENT MONITORS	LCD MONITORS
• Masimo	• Toshiba
• Spacelabs	• Sony
• GE	• Samsung
• Nihon Kohden	• Gateway
• Mindray	• Advantech
• Welch Allyn	• DELL
• Philips	• NEC
• Nonin	• LG

Options available for other manufacturers.
Please visit <http://www.amico.com/products/patient-monitor-mounts>

Product Classification



*RSS-PKG model shown here



*RSS-LCD model shown here

Package Model

RS_-42_-LCD/PKG (fill in blank with desired option)



5

No counterweight
max. 5 lbs (2.26 kg)



0

5" (127 mm) Shelf
w/ 30° tilt limit



2

6" (152 mm)
Drop on shelf



6

Base counterweight
max. 20 lbs (9.07 kg)



3

Philips IntelliVue Shelf
w/ 30° tilt limit



4

Accessory add
on top



7

Base and leg
counterweight
max. 25 lbs (11.33 kg)



6

5" (127 mm) Shelf for
back mounts



7

5" (127 mm) Shelf
w/slide-in VESA
75/100 Plate



8
3" (76 mm) Shelf
w/30° tilt limit



9
3" (76 mm) Shelf
for back mounts



W
Welch Allyn w/30°
tilt limit

*Weight ratings include 10 lbs in basket
*Highlight denotes counterweight placement

Standard Model

RS_-42_/LCD (fill in blank with desired option)



5

No counterweight
max. 10 lbs (4.53 kg)



0

5" (127 mm) Shelf
w/ 30° tilt limit



2

6" (152 mm)
Drop on shelf



6

Base counterweight
max. 25 lbs (11.33 kg)




3

Philips IntelliVue Shelf
w/ 30° tilt limit



4

Accessory add
on top



7

Base and leg
counterweight
max. 30 lbs (13.60 kg)



6

5" (127 mm) Shelf for
back mounts



7

5" (127 mm) Shelf
w/slide-in VESA
75/100 Plate



8
3" (76 mm) Shelf
w/30° tilt limit




9
3" (76 mm) Shelf
for back mounts



W
Welch Allyn w/30°
tilt limit

*Highlight denotes counterweight placement

Technical Drawings



Roll Stand Accessories



Tablet / Laptop Mount

- Adjusts to fit most tablet/laptop PCs

Please visit: www.amico.com/products/laptop-and-tablet-mounts



Cord Wrap

- Conveniently organizes cords

Part Number: CORDWRAP-R125



Soft Grip Handle

- Easy grab handle

Part Number: RS1-HANDLE



Power Bar

- Hospital grade, 4 power outlets
- 120 VAC, 60Hz

Part Number: PWRBAR4-R125



Power Holder

- Holds most power adapters for PCs, patient monitors and LCDs

Part Number: PWRHLD-R125



Wire Basket

- Holds up to 20 lbs (9 kg) (based on model)

Part Number: 100609-R125-BB



Device Integration Extenders

- Additional lightweight VESA mount

Please call to consult

Please visit www.amico.com for all technical specifications, installation instructions, detailed AutoCAD drawings and more Mounting Solutions.



Distributed By:

