


Philips Efficia & EarlyVue Monitor

Ordering Guide for Philips Sales



Device	Mounting Kits	Mounting Platform	Restrictions
 Philips EarlyVue VS30	Philips P/N: 989806108764 Amico P/N: PH-SSM-250-M-00 (Set) (Philips Efficia & EarlyVue Side to Side 9" Monitor Arm) Includes:  SSM-250 (Side to Side 9" Arm)  PH-EFF-PLT (Philips Efficia Mounting Plate)  MRS-COVER-019 (19" Mon. Channel Cover)  CRDHLDL (Dual Cord Holder)	 MRS Rail On Wall Mounted (sold separately) Amico P/N: MRS-S0(XX) (with different length) Philips P/N: 989806108759 19 – 19" Length Philips P/N: 989806108760 24 – 24" Length Philips P/N: 989806108761 34 – 34" Length Philips P/N: 989806108762 37 – 37" Length Philips P/N: 989806108763 48 – 48" Length	N/A

Pictures are not to scale.



Philips P/N: 989806108764

Amico P/N: PH-SSM-250-M-00 (Set)


Shipping Box: 20" (508 mm) x 12" (304.8 mm) x 6" (152.4 mm)

Weight: 6.5 lbs (2.95 kg) approx.

Philips Efficia & EarlyVue Monitor

Ordering Guide for Philips Sales



Device	Mounting Kits	Mounting Platform	Restrictions
 Philips Efficia CM10/12/100/120/150 Monitor CM10 Monitor (7.0 lbs [3.3 kg]) (7.6 lbs [3.6 kg]) with 3-cell battery (8.1 lbs [3.8 kg]) with 9-cell battery CM12 Monitor (11 lbs [5.0 kg]) (11.6 lbs [5.3 kg]) with 3-cell battery (12.1 lbs [5.5 kg]) with 9-cell battery CM100 Monitor (7.0 lbs [3.3 kg]) (8.1 lbs [3.8 kg])* with single battery CM120 Monitor (11.0 lbs [5.0 kg]) (12.1 lbs [5.5 kg])* with single battery (13.2 lbs [6.0 kg])* with dual batteries CM150 Monitor (14.5 lbs [6.7kg]) (15.6 lbs [7.2 kg])* with single battery (16.7 lbs [7.7 kg]) with dual batteries	Philips P/N: 989806108764 Amico P/N: PH-SSM-250-M-00 (Set) (Philips Efficia & EarlyVue Side to Side 9" Monitor Arm) Includes:  SSM-250 (Side to Side 9" Arm)  PH-EFF-PLT (Philips Efficia Mounting Plate)  MRS-COVER-019 (19" Mon. Channel Cover)  CRDHLDL (Dual Cord Holder)	 MRS Rail On Wall Mounted (sold separately) Amico P/N: MRS-S0(XX) (with different length) Philips P/N: 989806108759 19 – 19" Length Philips P/N: 989806108760 24 – 24" Length Philips P/N: 989806108761 34 – 34" Length Philips P/N: 989806108762 37 – 37" Length Philips P/N: 989806108763 48 – 48" Length	N/A

Pictures are not to scale.



Philips P/N: 989806108764

Amico P/N: PH-SSM-250-M-00 (Set)

Shipping Box: 20" (508 mm) x 12" (304.8 mm) x 6" (152.4 mm)

Weight: 6.5 lbs (2.95 kg) approx.