


# Philips IntelliVue

## Ordering Guide for Philips Sales



| Device                                                                                                                                                                                                                                                                                                                                            | Mounting Kits                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | Mounting Platform                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | Restrictions                                                                                                                            |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
|  <p><b>IntelliVue</b></p> <p><b>MX700/600</b> (18 lbs [8.2 kg])</p> <p><b>MX750</b> (22 lbs [10 kg])</p> <p><b>MX800</b> (26.4 lbs [12 kg])</p> <p><b>MX850</b> (24.3 lbs [11 kg])</p> <p><b>AD75</b> (22 lbs [10 kg])</p> <p><b>AD85</b> (24.3 lbs [11 kg])</p> | <p><b>Philips P/N: 989803224115</b><br/>Amico P/N: PH-AHM-053-M-01<br/>(Philips MX600-MX800, AD75 &amp; AD85 (Set) Adjustable Height Monitor Arm)<br/>Includes</p> <div>  <p><b>AHM-053 (PHILIPS)</b><br/>(Adjustable Height Arm)</p> </div> <div>  <p><b>MRS-COVER-019</b><br/>(19" Mon. Channel Cover)</p> </div> <div>  <p><b>AHM-HAND-00</b><br/>(AHM Handle)</p> </div> <div>  <p><b>DOWNPOST-09</b><br/>(Down Post 9")</p> </div> |  <p><b>MRS Rail</b><br/>On Wall Mounted<br/>(sold separately)</p> <p>Amico P/N: MRS-S0(XX)<br/>(with different length)</p> <p><b>Philips P/N: 989803224103</b><br/>19 – 19" Length</p> <p><b>Philips P/N: 989803224104</b><br/>24 – 24" Length</p> <p><b>Philips P/N: 989803224105</b><br/>34 – 34" Length</p> <p><b>Philips P/N: 989803224106</b><br/>37 – 37" Length</p> <p><b>Philips P/N: 989803224107</b><br/>48 – 48" Length</p> |  <p><b>Table Mount</b><br/>(supplied by Philips)</p> |

Pictures are not to scale.



**Philips P/N: 989803224115**

Amico P/N: PH-AHM-053-M-01 (Set)

Shipping Box: 24" (609 mm) x 11" (279.4 mm) x 6" (152.4 mm)

Weight: 13.75 lbs (6.24 kg) approx.