


Philips IntelliVue

Ordering Guide for Philips Sales



Device	Mounting Kits	Mounting Platform	Restrictions																		
 IntelliVue <table><tr><td>MX400</td><td>(13.3 lbs [6 kg])</td></tr><tr><td>MX450</td><td>(12.9 lbs [5.9 kg])</td></tr><tr><td>MX500</td><td>(12.6 lbs [5.7 kg])</td></tr><tr><td>MX550</td><td>(10 lbs [4.54 kg])</td></tr><tr><td>MP5</td><td>(13.45 lbs [6.1 kg])</td></tr><tr><td>MX700/600</td><td>(18 lbs [8.2 kg])</td></tr><tr><td>MX750</td><td>(22 lbs [10 kg])</td></tr><tr><td>MX800</td><td>(26.4 lb [12 kg])</td></tr><tr><td>MX850</td><td>(24.3 lbs [11 kg])</td></tr></table>	MX400	(13.3 lbs [6 kg])	MX450	(12.9 lbs [5.9 kg])	MX500	(12.6 lbs [5.7 kg])	MX550	(10 lbs [4.54 kg])	MP5	(13.45 lbs [6.1 kg])	MX700/600	(18 lbs [8.2 kg])	MX750	(22 lbs [10 kg])	MX800	(26.4 lb [12 kg])	MX850	(24.3 lbs [11 kg])	Philips P/N: 989806108781 Amico P/N: PH-SSM-053-M-00 (Set) (Philips MX400-MX850, MP5 Side to Side 12" Monitor Arm) Includes:  SSM-053 (Side to Side 12" Arm)   MRS-COVER-019 (19" Mon. Channel Cover) DOWNPOST-06 06" Down Post	 MRS Rail On Wall Mounted (sold separately) Amico P/N: MRS-S0(XX) (with different length) Philips P/N: 989806108759 19 – 19" Length Philips P/N: 989806108760 24 – 24" Length Philips P/N: 989806108761 34 – 34" Length Philips P/N: 989806108762 37 – 37" Length Philips P/N: 989806108763 48 – 48" Length	 Table Mount (supplied by Philips)
MX400	(13.3 lbs [6 kg])																				
MX450	(12.9 lbs [5.9 kg])																				
MX500	(12.6 lbs [5.7 kg])																				
MX550	(10 lbs [4.54 kg])																				
MP5	(13.45 lbs [6.1 kg])																				
MX700/600	(18 lbs [8.2 kg])																				
MX750	(22 lbs [10 kg])																				
MX800	(26.4 lb [12 kg])																				
MX850	(24.3 lbs [11 kg])																				

Pictures are not to scale.



Philips P/N: 989806108781

Amico P/N: PH-SSM-053-M-00 (Set)

Shipping Box: 20" (508 mm) x 12" (304.8 mm) x 6" (152.4 mm)

Weight: 6 lbs (2.72 kg) approx.