



| Device                                                                                                                                          | Mounting Kit                                                                                                                         | Mounting Platform                                                                                  | Restrictions | Add-ons                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|--------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  <p><b>Philips IntelliVue MX850</b><br/>24.3 lbs (11 kg)*</p> |  <p>N/A<br/>(Philips supplied Table Top Mounts)</p> |  <p>AHM-053</p> | N/A          |  <p>Philips IntelliVue X2</p>  <p>Multi-Measurement Module (MMS)</p>  <p>FMS-4 Module</p>  <p>FMS Module</p>  <p>Cable and Patient Leads</p> |

# Philips Ordering Guide Reference Chart



| Device                                                                                                                                                                             | Mounting Kit                                                                                                            | Mounting Platform                                                                                    | Restrictions                                   | Add-ons |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|------------------------------------------------|---------|
|  <p><b>Philips IntelliVue X2</b><br/>2 lbs (0.5 kg)*</p>                                          | <p>N/A<br/>(Mounts to the back of the MX850)</p>                                                                        | N/A                                                                                                  | N/A                                            | N/A     |
|  <p><b>Multi-Measurement Module (MMS)</b><br/>1.4 lbs (0.7 kg)*</p>                               | <p>N/A<br/>(Mounts to the back of the MX850)</p>                                                                        | N/A                                                                                                  | N/A                                            | N/A     |
|  <p><b>FMS-4 Module</b><br/>2.4 lbs (1.1 kg)*</p> <p><b>FMS Module</b><br/>7.7 lbs (3.5 kg)*</p> |  <p>N/A<br/>(Clamp from Philips)</p> |  <p>AHM-053</p> | <p>Clamps below with primary device on AHM</p> | N/A     |

# Philips Ordering Guide Reference Chart



| Device                                                                                                                  | Mounting Kit                                                                                           | Mounting Platform                                                                                  | Restrictions                            | Add-ons |
|-------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|-----------------------------------------|---------|
|  <p><b>Cable and Patient Leads</b></p> |  <p>CRDHLDR2</p>      |  <p>AHM-053</p> | Mounts below with primary device on AHM | N/A     |
|                                                                                                                         |  <p>101004-MON-BB</p> | N/A                                                                                                | Mounts on MRS rail with AHM             | N/A     |

\*Weight provided to help approximate mounting capacity. For reference only.