



## Features

- Allows healthcare professionals to customize the positioning of the sensor brackets according to the unique layout and requirements of each environment.
- This flexibility enhances usability and ensures optimal sensor placement for fall prevention monitoring.
- Fits Posey Sitter Elite Alarm 8208 wall bracket or 8276 bed bracket, or Medline MDT85WMT wall bracket.
- Can be mounted to any headwall, power column, flatwall, equipment boom, or direct rail option.
- Designed to meet infection control and prevention standards using non-porous, easy-to-clean materials (aluminum, steel, ABS plastic, nylon).
- Use CaviWipe or other disinfectant wipes (no abrasives, solvents, or trichloroethylene).
- Maximum load limit: 1.0 lb (0.45 kg).
- Warranty: 5 years


**NOTE:** The Posey Sitter Elite Alarm 8276 bed bracket is for illustration purposes only. Not included in the product.

## Key Dimensions & Weight

Part #	Adapter Option	Dimensions (W x D x H)	Weight	Image
PFS-PLT-A	Amico Horizontal Rail Mount (Compatible with Amico and Modular)	2.5" x 2.2" x 2.8" (6.4 cm x 5.6 cm x 7.1 cm)	0.6 lbs (0.3 kg)	
PFS-PLT-D	Dovetail Rail Mount (Compatible with Ohmeda Rail)	1.6" x 2.5" x 2.8" (4.1 cm x 6.4 cm x 7.1 cm)	0.5 lbs (0.2 kg)	
PFS-PLT-F	Fairfield Rail Mount (Push button optional)	1.5" x 0.8" x 2.9" (3.8 cm x 2.0 cm x 7.4 cm)	0.5 lbs (0.2 kg)	
PFS-PLT-G	Hill-Rom Vertical Rail Mount	1.5" x 2.6" x 2.8" (3.8 cm x 6.6 cm x 7.1 cm)	0.5 lbs (0.2 kg)	
PFS-PLT-MON	Monitor Rail Mount	2.5" x 0.5" x 4.5" (6.4 cm x 1.3 cm x 11.4 cm)	0.7 lbs (0.3 kg)	
PFS-PLT-PLH	DIN 25 mm (1") Rail Mount (Compatible with 0.75" [19.1 mm] to 1.5" [38.1 mm] horizontal pole)	1.6" x 2.4" x 6.2" (4.1 cm x 6.1 cm x 15.7 cm)	1.0 lb (0.5 kg)	

Part #	Adapter Option	Dimensions (W x D x H)	Weight	Image
PFS-PLT-PLV	With 0.75" (19.1 mm) to 1.5" (38.1 mm) Vertical Pole Mount	6.2" x 2.4" x 2.8" (15.7 cm x 6.1 cm x 7.1 cm)	1.0 lb (0.5 kg)	
PFS-PLT-R125	1 1/4" (31.75 mm) Pole Mount	2.8" x 1.8" x 2.8" (7.1 cm x 4.6 cm x 7.1 cm)	0.5 lbs (0.2 kg)	
PFS-PLT-R150	1 1/2" (38.1 mm) Pole Mount	3.1" x 2.2" x 2.8" (7.9 cm x 5.6 cm x 7.1 cm)	0.7 lbs (0.3 kg)	
PFS-PLT-R200	2" (50.8 mm) Pole Mount	3.5" x 2.6" x 2.8" (8.9 cm x 6.6 cm x 7.1 cm)	0.8 lbs (0.4 kg)	
PFS-PLT-S	HSI Mount-Vertical Rail Mount	1.5" x 2.9" x 2.8" (3.8 cm x 7.4 cm x 7.1 cm)	0.5 lbs (0.2 kg)	
PFS-PLT-V	Amico Vertical Rail Mount (Compatible with Modular Rail)	1.5" x 2.4" x 2.8" (3.8 cm x 6.1 cm x 7.1 cm)	0.5 lbs (0.2 kg)	

Please visit [www.amico.com](http://www.amico.com) for all technical specifications, installation instructions, detailed AutoCAD drawings, and more Mounting Solutions.



Distributed By:

