




Features

- 70" reach
- Ergonomic design for reduced clinician fatigue
- VESA 75/100 mounting pattern supports the wide range of devices sizing up to 24"
- Built-in cable management for hassle free operation
- Two extended arms provided to optimal reach and articulation
- Maximum Load: 2.21 to 14.33 lbs (1 to 6.5 kg)
- Designed to meet infection control and prevention standards using non-porous, easy-to-clean materials that are safe for use with hospital grade cleaning solutions

Specifications

Cable Routing

- USB-A, USB-C, PoE++ applications, Coaxial, Cat6, any cable with touch-safe voltage supply



WARNING: DO NOT ROUTE AC POWER CORD THROUGH THE ARMS.

Display Articulation

- Hinge Swivel: $\pm 185^\circ$
- Device Mount Rotation: $\pm 5^\circ$
- Viewing Angle Adjustment: -60° to $+20^\circ$

Display Device Support

- Mounting Pattern: 100 mm x 100 mm / 75 mm x 75 mm
- Mount: Single Device
- Display Type: Single monitor
- Panel Size Support: 15" to 24"
- Support Weight Range: 2.21 to 14.33 lbs (1 to 6.5 kg)

Height Adjustment

- Arm 1 Swivel (1): $\pm 90^\circ$
- Arm 1 Swivel (2): $\pm 185^\circ$
- Height Adjustment Range: 34.92" (88.7 cm)
- Tension Torque Adjust: Yes
- Arm Tilt (2): -55° to $+25^\circ$

Key Dimensions


Outside Dimensions (in extended position):

- Width: 7.87" (20.0 cm)
- Depth: 77.5" (196.8 cm)
- Height: 59.5" (151.2 cm)


Outside Dimensions (in folding position):

- Width: 7.87" (20.0 cm)
- Depth: 13.62" (34.6 cm)
- Height: 54.1" (137.4 cm)
- Physical Weight: 34.17 lbs (15.5 kg)

Figure 1: Long Reach Arm Articulation and Dimensions



Please visit **www.amico.com** for all technical specifications, installation instructions, detailed AutoCAD drawings, and more Mounting Solutions.



Distributed By: