

Vaccine Clinic

Mobile Supply & Administration Carts



Amico



Hummingbird **Laptop Workstation**

Non-Powered with Manual Lift

for Vaccine, Supply Storage and Administration

Highlights

- Facilitates easy access to vaccine administration supplies
- Secure laptop compartment and optional integrated RFID badge reader
- Easily connect and mount accessories, such as barcode scanners and sharps containers
- Durable surfaces suitable for cleaning with Virox, CaviWipes, and other common cleaners between patients
- Easily remove accessories for terminal cleaning and sterilization
- Future-Proof – ability to add on accessories and modify cart in the future
- Optional on-board power (6,12, or 18 hour rated battery)
- Optional LCD/AIO mount with CPU compartment

As Shown

- Non-Powered, manual lift (with pneumatic assistance)
- Slide-out Keyboard Tray
- Sharps Container Holder
- Barcode Scanner Dock Mount
- Glove Box Holder
- Sanitizer Pump Mount
- 4-drawer, Single Bank Storage



Accessories (See Page 4)



Hummingbird Lite **Laptop Workstation**

Non-Powered with Manual Lift

for Patient Screening, Registration, and Administration

Highlights


- Lightweight for fast and easy deployment throughout the hospital and remote clinic
- Lockable laptop security bracket ensures assets stay with the cart
- Durable surfaces suitable for cleaning with Virox, CaviWipes, and other common cleaners
- Easily remove accessories for terminal cleaning and sterilization
- Future-Proof – ability to add on accessories and modify cart in the future

As Shown

- Laptop model with Integrated Security Bracket (laptop not included)

Options

- Laptop, Tablet, or LCD/AIO
- Integrated Security Bracket



RST Cart

Tablet Roll Stand

for Patient Screening and Registration

Highlights

- Lightweight, secure solution
- Tablet adjustment: tilt, pan, and rotate for optimal ergonomic positioning and ease of use
- Durable surfaces suitable for cleaning with Virox, CaviWipes, and other common cleaners
- Easily remove accessories for terminal cleaning and sterilization
- Future-Proof – ability to add on accessories and modify cart in the future

As Shown

- 42" (106.7 cm) tall with LCD Mounting Plate to attach Secure Tablet Enclosure
- Tablet Enclosure
- Handle
- Sturdy 24" (53.3 cm) Base with Locking Casters
- 10" x 6" x 9" (25.4 x 15.2 x 22.9 cm) Wire Basket with Cordwrap

Options

- Gooseneck Extension with Locking Tablet Holder
- 5" x 5" x 4" (12.7 cm x 12.7 cm x 10.2 cm) Wire Basket with Cordwrap



Accessories for Workstations



Glove Box Holder

Storage for boxes of gloves
(works with tissue boxes too)



Hand Sanitizer Holder

Securely store your hand sanitizer pump




Wipes Holder

Plate for mounting your sanitizer dispenser



Cordwrap

Neatly store your cables and cords when not in use. Attaches to basket or cart.



Storage Baskets

Sterilizable wire baskets for storing items such as masks, gowns, and supplies
Small 5" x 5" x 4" (12.7 cm x 12.7 cm x 10.16 cm)
Medium 10" x 6" x 9" (25.4 cm x 15.24 cm x 22.86 cm)



Drawers

Storage for vaccines and supplies 4x slim or 1x wide drawers up to 3 banks



RFID

RFID Integration facilitates secure PC access and enhanced infection control



Printer Mount

Mounting plate accommodates commonly used small printers including Zebra QL220



Shelf

Shelf for placing accessories such as printers
Small 6" x 6" (15.24 cm x 15.24 cm)
Medium 12" x 12" (30.48 cm x 30.48 cm)
Large 12" x 20" (30.48 cm x 50.80 cm)



Barcode Scanner

Mount your wired or wireless barcode scanner beside your display or worktop



Universal Sharps Container

Adjustable mount to accommodate common sharps containers

Amico Accessories Inc. | 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada
Toll Free Tel: 1.877.264.2697 | Tel: 905.763.7778 | Fax: 905.763.8587
Email: accessories@amico.com | www.amico.com

Represented By: