




Wall-Mounted Computer Workstations



Amico



At Amico, we believe in improving healthcare for every stakeholder. The heroes in our story are Patients, Caregivers, Technicians, and Management – *they* are the leaders of our development process and the reason you find Amico innovations throughout the healthcare environment. *Our* expertise is in listening to them.





Patients

Elevated Healthcare

Support patient-centered care and improved communication at the bedside. Optimized patient experiences increase compliance to treatment plans, reduce return visits, and most of all - improve outcomes.



Improved
Communication



Improved Experience



Increased Treatment
Compliance



Better Outcomes



Caregivers

Effortless Workflow

Designed collaboratively with caregivers to solve *real* workflow challenges across a variety of clinical environments. Ergonomics, ease-of-use, and durability are combined to achieve functional workstations that facilitate focus on patient care.



Increased Workflow
Efficiency



Improved Charting and Medication
Administration Accuracy



Less Strain



More Energy



Technicians

Peace of Mind

Durable, infection-control friendly designs with integrated cable management, including thoughtful features to maximize uptime, and minimize installation and service requirements. A suite of training and installation options available to suit the needs of your facility.



Ships Assembled



Minimal Maintenance



Easy to Clean



Multi-level Support



Management

Productivity and Savings

Efficiency, durability, and adaptability translate into long-term productivity and savings across the facility. Amico computer workstations and accessories are a clinically optimized solution that minimizes installation and service with easy integration into our full suite of products and equipment.



Installation Savings



Better Patient Outcomes



Increased Caregiver Productivity



Part of Amico Ecosystem


Falcon

LCD and Keyboard Workstation





Anatomy of the Falcon



Falcon Models

Without SideKick



With SideKick



With SideKick and Extension
6" (15.2 cm) - 16" (40.6 cm) available



SideKick Secondary Worksurface



Add optional RFID Integration for secure access to your workstation



SideKick Secondary Worksurface gives caregivers room to work, an integrated light for night charting, and optional RFID security



Flip up your SideKick Secondary Worksurface for quick access to storage and magnetic whiteboard


Condor

Adjustable Height Channel Workstation





Anatomy of the Condor



Condor Models

Ultra Flush



Flush with SideKick



With SideKick and Extensions
9" (22.9 cm) - 24" (61 cm) available



SideKick Secondary Worksurface



Add optional RFID Integration for secure access to your workstation



Flip up your SideKick Secondary Worksurface for quick access to storage and magnetic whiteboard



Condor Lite available with fixed height

Eagle

Dual Arm Workstation



Heron

Tablet Mount



Hawk

Laptop Mount



Raven

Overhead-Mounted LCD &
Keyboard Workstation



Sparrow Monitor Arm

Light-Duty



Anesthesia Cart Mounts



Workflow

Mounting Locations

Amico Computer Workstations and Accessories can be mounted in a range of locations and to a variety of equipment rails and headwalls. Order your Computer Workstation and Accessories with a mounting rail, or a mounting adaptor for seamless integration into your existing architecture.



Ergonomics

Ease of use and ergonomic adjustability are fundamental design aspects of our products. Each product is developed through our healthcare facility-lead design process resulting in user-friendly, workflow-friendly equipment.



Workflow

Installation & Maintenance

Low maintenance is consistently a high priority with our healthcare partners. We've designed our products with minimal installation preparation, zero assembly, and easy-to-clean components, reducing workload on the facilities team.

Remote and on-site training, as well as white glove installation services available.



Amico Ecosystem

Amico Wall-Mounted Workstations, Arms, Accessories, and many of our other products are made to be moved from room to room, enabling low cost adjustments for changes in acuity and workflow.



Product Features

	Falcon <i>Combo Arm</i>	Falcon Lite <i>Side-to-Side Combo Arm</i>	Condor <i>Adjustable Height Channel</i>	Condor Lite <i>Fixed Height Channel</i>	Eagle <i>Dual Arm</i>	Raven <i>Combo Arm</i>	Hawk <i>Single Arm Laptop Mount</i>	Heron <i>Tablet Mount</i>	Sparrow Light- Duty Arm <i>Monitor Arm</i>
Features	● Standard				○ Optional				
Adjustable Height	●		●		○	●	○	○	●
Display Tilt, Pan, Rotate Adjustments	●	●	●	●	●	●		○	●
Independent Display- Height Adjustment	●	●			●	○			
Side-to-Side Movement	●	●	●	●	●	●	●	●	●
SideKick Secondary Workspace	○	○	○	○	○	○			
Add-on Extension Optional for All	●	●	●	●	●	●	●	●	○
Easy Relocation to Other Mounting Locations	●	●			●	●	●	●	
Integrated Nightlight	○	○	○	○	○	○	○	○	
Integrated RFID	○	○	○	○	○	○			
Integrated Cable Management	●	●	●	●	●	●	●	●	●
Compact Folding	●	●	●	●	●	●		○	
Add-on Accessories	○	○	○	○	○	○	○	○	○
Color: White	●	●	●	●	●	●	●	●	
Color: Black	○	○	○	○	○	○			●
Hallway-friendly Rating	Better	Better	Best	Best	Good	Good	Good	Good	

Customize Your Workflow

Options



Keyboard Trays

Magnetic keyboard tray also available as a factory installed add-on – part number: KEYBRDD-MAGPLT




Monitor Arm Colors



Channel Covers



Writing Surfaces



SideKick
Secondary Worksurface




Dual Display



CPU Mount - Large
Mount for larger form-factor PC's
behind your display



CPU Mount - Small
Mount for small form-factor PC's
behind your display



Thin Client Mount
Behind LCD




Enclosed CPU Mount



Tangent Rugged Mini Mount



Extension Options
Various extension sizes available



Articulating Shelf
Shelf on articulating arm can
move side-to-side



Universal Scanner Holder



Barcode Scanner Mount
Horizontal mounting plate for
your barcode scanner



Power Supply Mounts
Holds most power adapters for
laptops, patient monitors and LCDs



RFID Reader
Dual-frequency versatile RFID
reader, easy to install and clean

See spec sheets for option availability by product

Customize Your Workflow

Accessories



Glove Box Holder

Easy to load and mounts to any accessory bar or wall channel



Hand Sanitizer Holder

Ideal for personal sized containers



Universal Sharps Holder

Ideal for mini sharps box and universally compatible



Hospital Grade Power Bar

Hospital grade, 4 power outlets






Chart Holder

Locate your charts close at hand, beneath the worktop



Screen Handles

Ergonomic handles make adjusting your display easy



Antimicrobial Keyboard and Mouse

(Ask sales rep for options)




Wrist Rest

Antimicrobial comfort wrist rest and Virox friendly



Mouse Pad

Wipeable self adhering and non-reflecting, ideal for optical mice



Mouse Holder

Neatly stores mouse when not in use



Baskets

Sterilizable wire baskets
10" x 6" x 9" (25.4 cm x 15.2 cm x 22.9 cm)
10" x 10" x 4" (25.4 cm x 25.4 cm x 10.2 cm)
5" x 5" x 4" (12.7 cm x 12.7 cm x 10.2 cm)



Cuff Baskets

Sterilizable wire baskets
13" x 9" x 7" (33.02 cm x 22.9 cm x 17.8 cm)



Need a **Mobile** Workstation?
Ask us about the
Hummingbird



Amico Accessories Inc. | 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada
Toll Free Tel: 1.877.264.2697 | Tel: 905.763.7778 | Fax: 905.763.8587
Email: accessories@amico.com | www.amico.com

Represented By: