

ICU Series B-AM1-400802 Bed



Key Features:

- 550 lbs (250 kg) weight capacity
- Standard 5th wheel for improved mobility
- 21" (53 cm) low height
- Storage area in base for two oxygen E-cylinder tanks and belongings
- Radiolucent fowler section
- Dual brake pedal mechanisms
- Adjustable knee gatch for enhanced comfort
- Fold-down siderails
- 8" (20 cm) wheel casters
- Pressure relief surfaces

Standard Features:

- LINAK™ DC Actuators
- Modular battery backup
- LED night light
- Four IV pole sockets
- Four integrated roller bumpers
- Metal litter deck
- X-ray cassette holder
- CPR release system
- Trendelenburg and reverse Trendelenburg
- Patient controls on inner siderails
- Flat litter deck for compatibility with standard patient surfaces
- Adjustable bed height
- Adjustable foot extension
- Adjustable litter with chair position
- Integrated digital scale
- 3-Zone Bed exit alarm
- Scale and fowler angle indicator on siderail
- Foley bag holder
- Integrated push handles
- 5th wheel for improved mobility
- Brake not set alarm

The new Amico ICU bed features a versatile design suitable for a wide range of patients. The ICU bed is differentiated by its enhanced chair position as well as its ergonomic and robust design. This product truly helps provide a more meaningful level of care to patients.

This bed has a battery back-up system to enable motion even when the bed is disconnected or during patient transport. The ICU bed also features LINAK™ DC actuators that offer a reliable yet silent motion for the bed's electromechanical system.

Technical Specifications



Model Number

Model Number	B-AM1-400802
--------------	--------------

Overall Dimensions

Overall Length	
Retracted	83" (211 cm)
Extended	91" (231 cm)
Overall Width	
Siderails Up	42" (107 cm)
Siderails Down	38" (96.5 cm)
Patient Surface	
Retracted	32" x 71" (81 cm x 180 cm)
Extended	32" x 79" (81 cm x 201 cm)
Height Range (to the surface of the litter)	
High	38" (96.5 cm)
Low	20.5" (52 cm)

Litter Positioning

Trendelenburg/Reverse Trendelenburg	± 18°
Backrest	0° - 70°
Knee Gatch	0° - 21°

Weight Capacity & Caster Diameter

Weight Capacity	550 lbs (250 kg)
Caster Diameter	6" (15 cm)

Electronics

Volts	
Domestic	120 VAC
International	220 VAC
Ampere Rating	
Domestic	7.5 A
International	4.0 A
Current Leakage	<300 µA
Frequency	50 - 60 Hz
Hospital Grade Plug	Yes

Optional Features & Accessories:

- Pressure relief surfaces
- Removable IV pole - 2 stage
- IV hook
- Siderail covers
- Overbed table
- Communication to nurse call systems
- Oxygen tank holder
- Hand pendant controls

Amico Beds Corporation reserves the right to change product specifications without notice.

Warranty: One year parts only.

***Do not turn the bed on its side to move through doorways. Doing so may cause loss of functionality or damages to the bed.**