


M-THEA-WL003-BXXXXX

Non-Powered Mattress Behavioural



Key Features:

- High density foam with superior envelopment provides comfort to the patient
- Softer supportive comfort layer
- Enhanced low slip bottom prevents shifting of the mattress when it is placed on the bed frame
- Lightweight design allows for easy caregiver functionality
- Inside seams, welded closed - no zippers
- 4-way stretcher cover that is fluid resistant, fire retardant, antimicrobial and antifungal
- Fire Barrier Sleeve
- Part number: M-THEA-WL003-BXXX

Amico Beds Corporation reserves the right to change product specifications without notice.

Amico Beds Corporation | 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada | 600 Prime Place, Hauppauge, NY 11788, USA
Toll Free Phone: 1.877.462.6426 | Toll Free Fax: 1.866.440.4986 | Tel: 905.764.0800 | Fax: 905.764.0862 | www.amico.com

AB-SP-M-THEA-WL003-BXXXXX-NON-POWERED-MATTRESS-BEHAVIOURAL 04.14.2026