



Key Features:

- Integrated scale system
- Centrally controlled brake system
- Collapsible siderails
- Dual side foot pedals
- Retractable fifth wheel steering system
- X-Ray compatible backrest
- Shock absorbing, non-marking bumpers
- Foley bag holder
- Ergonomic push handles
- Storage area in base for oxygen E-cylinder tank
- Permanent IV pole

The Amico Hydraulic Patient Transfer Scale Stretcher (S-AM1-3006) provides caregivers the control they need in the emergency environment. The stretcher features an integrated scale system, allowing caregivers to weigh patients. Additionally, the patient exit alarm function detects when a patient exits the stretcher and gives an audible alarm when that happens. The stretchers can be adjusted to various positions – raising, lowering, Trendelenburg, and reverse Trendelenburg – by actuating the hydraulic pistons. The back section is adjustable through the use of gas springs and the knee gatch can be raised and lowered with a manual crank.

The stretcher is made of carbon steel and painted by an electrostatic spray (powder-coat). The base cover is made from injection molded ABS plastic and the locking siderails are made from a combination of steel, aluminum extrusion, and plastic parts for increased durability and to facilitate easy cleaning. The back section includes an X-Ray compatible board.

The stretcher comes complete with an IV pole, bumper wheels, and actuator dust covers to protect the hydraulic actuators.

Technical Specifications



Model Number

Model Number	S-AM1-3006
--------------	------------

Overall Dimensions

Overall Length	84.5" (214 cm)
Overall Width	33" (84 cm)

Litter Positioning

Trendelenburg/Reverse Trendelenburg	16°
Maximum Fowler Section Inclination	80° ± 5
Maximum Knee Section Inclination	35° ± 5

Height Range

Lowest	22" (56 cm)
Highest	34" (86 cm)

Surface Dimensions

Surface Width	25" (64 cm)
Surface Length	75" (191 cm)
Surface Thickness	4" (10 cm)

Electronics

Volts and Amps	
Domestic	120 VAC - 4A
International	220 VAC - 2.9A
Frequency	50 - 60 Hz

Weight Capacity & Caster Diameter

Maximum Weight Limit	550 lbs (250 kg)
Caster Diameter	8" (20 cm)

Optional Features & Accessories:

- Utility tray
- Overbed tables
- 4" (10 cm) mattress
- Four sided brake
- Removable IV pole

Amico Beds Corporation reserves the right to change product specifications without notice.

Warranty: One year parts only.

