

# Titan Series OB/GYN Stretcher



## Key Features:

- 550 lbs (250 kg) weight capacity
- 24" (61 cm) low height
- Storage area for patient belongings
- Easily accessible brake pedals
- Fold-down siderails
- 8" (20 cm) wheel casters
- Pressure relief surfaces
- Easy to operate lowering leg section
- Retractable fluid basin holder

## Standard Features:

- Two IV pole holders
- Roller bumpers
- Adjustable fowler section
- Fifth wheel steering
- Two side mounted pedals for patient surface height adjustment
- Hydraulic lift system capable of supporting up to 550 lbs (250 kg)
- Central brake locking system
- Ergonomic push handles
- Side bumpers
- Head and foot-end mattress retainer
- Retractable fluid basin holder
- Adjustable calf supports
- Pull-out stirrups

The Amico Titan OB/GYN Stretcher is designed with functionality and convenience in mind, allowing caregivers to provide a more meaningful level of care to patients.

This product offers key features such as a one hand operated lowering leg section, a deep well storage area for patient belongings and an ergonomic design to help improve efficiency of healthcare. To accommodate obstetrics patients with limited mobility, the Titan OB/GYN Stretcher is built on Amico's best-in-class, low height stretcher platform.

The Titan OB/GYN Stretcher can also serve as an interim transport stretcher for moving patients between departments.

# Technical Specifications



## Model Number & Approval

Model Number	S-AM1-3005
Agency Approval	cULus

## Overall Dimensions

Overall Length	77" (195.5 cm)
Overall Width	38" (96.5 cm)

## Litter Positioning

Trendelenburg/Reverse Trendelenburg	16°
Maximum Fowler Section Inclination	80°

## Height Range

Lowest	24" (61 cm)
Highest	33" (84 cm)

## Surface Dimensions

Surface Width	30" (76 cm)
Surface Length	71" (180 cm)
Surface Thickness	4" (10 cm)

## Weight Capacity & Caster Diameter

Weight Capacity	550 lbs (250 kg)
Caster Diameter	8" (20 cm)

## Optional Features & Accessories:

- 4"/5" (10 cm/12.7 cm) high density foam mattress with or without fire barrier
- Foley bag holder
- Overbed tables
- Permanent IV pole
- Calf supports
- Pull-out stirrups

Amico Beds Corporation reserves the right to change product specifications without notice.

**Warranty:** One year parts only.