

Titan Series Transport Stretcher



Key Features:

- 550 lbs (250 kg) weight capacity
- Standard fifth wheel for improved mobility
- 24.25" (61.5 cm) low height
- Storage area in base for two oxygen E-cylinder tanks and belongings
- Radiolucent fowler section
- Dual brake pedal mechanisms
- Adjustable knee gatch for enhanced comfort
- Fold-down siderails
- 8" (20 cm) wheel casters
- Pressure relief surfaces

Standard Features:

- Four IV pole holders
- Roller bumpers
- Adjustable fowler section
- Fifth wheel steering
- Fold-down siderails
- 8" (20 cm) wheel casters
- Hydraulic lift system
- Central brake locking system
- Ergonomic push handles
- Side bumpers
- Head and foot-end mattress retainer
- Integrated foot push handle

The Amico Titan Transport Stretcher uses an innovative design and incorporates distinctive features aimed to help caregivers provide a more meaningful level of care to patients.

This product is specifically designed to meet the needs of customers by offering key characteristics such as a dual brake system, deep well storage area in the base and an ergonomic design to help improve efficiency of healthcare. With the heightened awareness regarding fall prevention, Amico is proud to offer a Transport Stretcher that is the best-in-class low height stretcher in the market.

Technical Specifications



Model Number

Model Number	S-AM1-3000
Agency Approval	cULus

Overall Dimensions

Overall Length	84.5" (214 cm)
Overall Width	33" (84 cm)

Litter Positioning

Trendelenburg/Reverse Trendelenburg	14°
Maximum Fowler Section Inclination	75°
Maximum Knee Section Inclination	35° ± 5

Height Range

Lowest	24.25" (61.5 cm)
Highest	34.5" (88 cm)

Surface Dimensions

Surface Width	25" (63.5 cm)
Surface Length	75" (191 cm)
Surface Thickness	4" (10 cm)

Weight Capacity & Caster Dimension

Maximum Weight Limit	550 lbs (250 kg)
Caster Diameter	8" (20 cm)

Optional Features & Accessories:

- 4"/5" (10 cm/12.7 cm) high density foam mattress with or without fire barrier
- Utility tray
- Foley bag holder
- Overbed tables
- Permanent IV pole

Amico Beds Corporation reserves the right to change product specifications without notice.

Warranty: One year parts only.