



MAJESTIC SERIES HORIZONTAL HEADWALL SYSTEM

DRAWING # 0X

TWO TIER WITH NO CHASE
(M/N: HW00-02T-P-M084)



CEILING LINE



FINISHED FLOOR LINE

VIEW A-A

IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT / SIGNATURE


DATE

PHONE NO.

IF MIRRORED UNITS ARE REQUIRED, THE CUSTOMER MUST INDICATE THE QUANTITY OF UNITS BELOW AND AMICO WILL PROVIDE A SEPARATE DRAWING SHOWING THE MIRRORED LAY OUT AND QUANTITY QUANTITY MIRRORED: _____

TYPE:
QUANTITY:

SYSTEM DETAILS		
SYMBOL	QTY.	DESCRIPTION
HC	4	HORIZONTAL EQUIP. TRACK
O	2	AMICO GAS. OXYGEN
A	2	AMICO GAS. MED AIR
V	4	AMICO GAS. VACUUM
E2	8	RECEPTACLE - DUPLEX RED
N2	2	RECEPTACLE - DUPLEX IVORY
NC	1	PROVISION - NURSE CALL 3G
CB	1	PROVISION - CODE BLUE 1G
DA	2	PROVISION - DATA 1G
SLM	4	VACUUM SLIDE, TRACK MTD.



Amico
85 Fulton Way
Richmond Hill, Ontario
L4B 2M4, CANADA
Tel: (905) 1-866-440-4986
Fax: (905) 764-0802
www.amico.com

HOSPITAL	HOSPITAL
LOCATION	LOCATION
QTY.	UNITS AS SHOWN

A. NURSE CALL MFGR.:	MODEL NO.:	DRWG. NO.:
B. MEDICAL GAS MFGR.:	TYPE CONNECTION:	YYMMDD-XXXX
C. FINISH:	CEILING HEIGHT:	DRAWN BY: XX
		CHECKED BY: XX
		REV. NO.: 00
		DATE: MM/DD/YY