

CEILING COLUMN OUTLET DISS

MODEL NUMBERS:

The XXX defines the Gas:

Oxygen	= OXY
Medical Air	= AIR
MedVac	= VAC
Nitrous Oxide	= N2O
Nitrogen	= NIT
Carbon Dioxide	= CO2
Carbogen	= CAR
Waste Anesthetic Gas Disposal	= WAG (NFPA)






O-DISCCOL-L-XXX

The L defines the Language:

U = English (NFPA) Oxy, Air, Vac and WAG only
S = Spanish (NFPA) Oxy, Air, Vac and WAG only

Example:

English NFPA - DISS, Oxygen = O-DISCCOL-U-OXY




85 Fulton Way
Richmond Hill, Ontario
L4B 2N4, CANADA

Tel: (905) 764-0800
Fax: (905) 764-0862

www.amico.com