


# WALL OUTLET QUICK DISCONNECT PURITAN - BENNETT COMPATIBLE



1/2" O.D. (12.7mm) (3/8" nominal) size  
type "K" copper pipe, rotates 360 deg.  
for entry from any angle



### MODEL NUMBERS:

The XXX defines the Gas:  
 Oxygen =OXY  
 Nitrous Oxide =N2O  
 MedVac =VAC  
 Waste Anesthetic Gas Disposal =WAG

O-PBWAL-L-XXX

The L defines the Language:  
 U =English (USA Colors)  
 Oxy, Air, Vac, and WAG only  
 S =Spanish (USA Colors)  
 Oxy, Air, Vac, and WAG only

Example:  
 OXYGEN NFPA = O-PBWAL-U-OXY



Note:  
 Amico recommends that you have a blank space  
 between outlets and slides.  
 Maximum 3 outlets may be ganged together  
 without additional support.

85 Fulton Way  
 Richmond Hill, Ontario  
 L4B 2N4, CANADA  
 Tel: (905) 764-0800  
 Fax: (905) 764-0862  
[www.amico.com](http://www.amico.com)