

FlexCam 4K

Specifications



Description

The FlexCam 4K Camera is designed to provide clinicians with superior visuals to support big screen surgeries. This 4K Camera achieves such graphics through its exceptional functions which include focus, zoom, brightness, rotation, freeze and source select. Mounting directly to the Amico Spring Arm, the 4K Camera does not require a light head and can complete a 360-degree rotation with ease.



Built-in Control Pad



Removable Sterilizable
Camera Cover

Amico's FlexCam 4K Camera yoke can be controlled both manually and motorized for adjustments.