


Operating Room Solutions



Amico

Amico Group of Companies

Amico Corporation, founded in 1974, designs, manufactures, and markets a broad range of products for virtually every department in the health care facility. These products are manufactured in eight facilities in the U.S. and Canada, which occupy over **800,000** square feet.



Amico Group of Companies - One Point of Contact. One Total Solution.

Amico in the OR

Amico Clinical Solutions creates state-of-the-art products to ensure the greatest possible clinical experience. Amico is committed to providing customized solutions tailored to each surgeon's needs. Through close collaboration, Amico can design and deliver solutions that seamlessly integrate into your existing or new OR space. Amico's advanced technology optimizes OR efficiency and ensures the best possible OR outcomes.

Your Total OR Solution



Whether a new build or a renovation project, Amico Clinical Solutions can help customize your OR. Amico's innovative solutions can transform your facility, and ensure optimal clinician and patient OR experience.

- | | | |
|------------------------------|------------------------------|----------------------------------|
| ① Equipment Transport System | ⑧ Spring Arm Monitor Holders | • Video Integration (not shown) |
| ② GoLift | ⑨ Equipment Boom | • Pipeline Equipment (not shown) |
| ③ Patient Lift Pendant (PLP) | ⑩ Flowmeters | - Alarms |
| ④ Anesthesia Boom | ⑪ Medical Gas Outlets | - Medical Air Systems |
| ⑤ iCE Surgical Lights | ⑫ Vacuum Regulator | - Medical Vacuum Systems |
| ⑥ 4K Camera | | - Medical Gas Manifold |
| ⑦ Large Display Monitor | | - Valves |

iCE30m Surgical Lighting System



Shadow Reduction

Offers precise focus and optimal illumination to the surgical area, effectively eliminating shadows and ensuring maximum performance.



Sterile Field Control

Capacitive touch handle to allow for sterile field control. Manual focus without moving the lighthead - adjust the focal point up or down with one twist of the handle.



Multi-Color

Ranging from warmest to coldest white, surgeons may switch between four temperatures. Color temperatures include 3,500 K, 4,000 K, 4,500 K, and 5,000 K.



Ease of Movement

Dual yoke system to allow for flexible and easy positioning of the light.

Video Integration / Cameras



Integrated HD Camera

All cameras offer functions such as focus, zoom, brightness, rotation, freeze, and source select.



Integrated 4K Camera

Option of 4K wired or 4K wireless camera, both provide integrated functions such as optical 20x zoom, digital 12x zoom, and auto focus for a superior image.



Spring Arm Camera

Provides clinicians with superior visuals to support big screen surgeries.



Supports Big Screen Surgery

Large display monitor allows for vertical and horizontal movement with multiple height adjustment options.

iCE OR Booms



Anesthesia Boom

Equipped with fully customizable features such as medical gas outlets, electric outlets, monitor arms, and integrated accessory rails.



Equipment Boom

Customizable console and shelves allowing for integration of multiple accessories.



Hammerhead Boom

Its horizontal design is perfect for tight spaces and provides health care professionals with more space to work under the console or over top of the anesthesia cart, reducing its footprint.



Concealed Hinged Back Door

Allows for easy maintenance of the boom and future proof capability.

Equipment Transport System (ETS)



Streamlined Patient Transport

Hands-free transport over uneven floors and elevator thresholds allowing for safe transfer from OR to ICU.



Boom to Roll Stand

Easily lift IV pole from roll stand to boom with effortless foot pedal without relying on battery power.



Stacking Roll Stand

Stackable roll stands offer a space-saving solution.



Boom to ICU Bed

IV rack can be directly mounted to the bed from the boom without need for the roll stand.

ice 30^m LED OR Configurations

With so many configuration options to choose from, you can customize your Operating Room to maximize space and create a more efficient work environment.




SERIES 1 Single mount, single light



SERIES 1 Single mount, dual lights



SERIES 1 Single mount, 32" monitor



SERIES 1 Single mount, single light, 32" monitor



SERIES 3 Triple: Dual lights, 32" monitor on satellite



SERIES 3 Quad: Dual lights 32" monitor, 32" monitor on satellite



Amico Clinical Solutions Corp. | 122-B East Beaver Creek Road, Richmond Hill, ON L4B 1G6, Canada
600 Prime Place, Hauppauge, NY 11788, USA
Toll Free Tel: 1.877.462.6426 | Tel: 905.764.0800 | Fax: 905.764.0862
Email: info@amico.com | www.amico.com