

GoLift

Specifications



Product Features:

- Wired Hand Control
- Quick Connect Carry Bar
- Magnetic Hand Control
- Light Weight and Portable
- LED Indicators
- Infection Control
- Durable and Strong
- Easy to Service

Safety Features:

- Emergency Stop
- Emergency Lowering Device
- Upper Limit Detection
- Lower Limit Detection
- Slack Tape Sensor
- Free Fall Brake (over speed governor)
- Low Battery and Dead Battery Alarms
- Soft Start and Stop
- Overload Protection
- Emergency Manual Lowering

Single Lift GO-LFT-400

GoLift 400
Lifts up to 400 lbs / 181 kg



Single Lift GO-LFT-700

GoLift 700
Lifts up to 700 lbs / 318 kg



Tandem Lift GO-LFT-1000

GoLift 1000
Lifts up to 1000 lbs / 454 kg



Specifications:

Maximum Lifting Speed

No load:	2 inches/second
GoLift400:	1 inch/second
GoLift700:	1 inch/second
GoLift1000:	1 inch/second

Unit Weight

GoLift400:	10 lbs (4.5 kg)
GoLift700:	10 lbs (4.5 kg)
GoLift1000:	20 lbs (9 kg) <i>(Tandem Lifts)</i>

Safe Working Load (SWL)

GoLift400:	400 lbs (181 kg)
GoLift700:	700 lbs (318 kg)
GoLift1000:	1000 lbs (454 kg)

Dimensions

GoLift400 / GoLift700:	
• Length:	7.75" (196.8 mm)
• Width:	7.75" (196.8 mm)
• Height:	4.25" (107.9 mm)
• Strap Length:	84" (2133.6 mm)
GoLift1000:	
• Length:	15.5" (393.7 mm)
• Width:	7.75" (196.8 mm)
• Height:	4.25" (107.9 mm)
• Strap Length:	84" (2133.6 mm)

Number of Lifts Per Charge (Duty: 10/90)

25% of Strap at Midrange ¹ :	
• 375 with 185 lbs (84 kg)	
• 170 with 400 lbs (181 kg)	
• 102 with 700 lbs (318 kg)	
• 140 with 1000 lbs (454 kg)	
Charging time:	2 - 4 hours
¹ Calculated using 5Ah Packs	

Lift Case

Flame retardant ABS

Batteries

High Capacity, Nickel Metal Hydride (NiMH)	
Standard:	2x 14.4V (3.3Ah)
Optional:	2x 14.4V (5Ah)

Battery Charger

Input:	100-240V AC
Output:	36V DC, 1.0A, 40W

Hand Control

Protection Class: IPX4

Maximum Sound Level

<55 dB


Expected Service Life

10 years or 22,500 cycles


Approvals

Certified to:	CAN/CSA-C22.2 No 60601-1:08 UL 60601-1:08
Tested to:	ISO 10535-06, CE


GoTrack Configurations:



Regular



Double



Triple

