

# GoLift Scale

## Specifications



### General Specifications:

The scale provides an integrated method to the patient lift system to measure a patient's weight. It is a modular component that is mounted between the bottom of the strap and the carry bar. It features the Amico quick release system so that it can be removed and attached without the use of tools. Since the scale can be easily removed, it can be moved around from room to room as required.

### Specifications:

Two versions available:

- Class II - Displays units in kg or lbs, 1000 lbs capacity
- Class III - Displays units in kg, 400 kg capacity, 2014/31/EU - CE conformity

### Features:

- Over capacity sensor
- Weight positive lock (to prevent accidental detachment)
- Runs off a 9V battery with low battery indicator
- Smart stable weight function, battery saving function, and zeroing option.
- Provides up to 3000 weight readings before needing replacement
- Sealed casing for infection control (clean using mild soap and a soft damp cloth with water).