


# Sling Details and Sizing

## Guidelines For Choosing Your Sling Sizes



|                                     | S                                     | M                                     | L                                     | XL                                    |
|-------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|
| <b>GoBasic</b><br>With Head Support | Height 29"<br>Shoulder 20"<br>Hip 16" | Height 30"<br>Shoulder 24"<br>Hip 20" | Height 31"<br>Shoulder 28"<br>Hip 24" | Height 32"<br>Shoulder 32"<br>Hip 28" |
| <b>GoBasic</b><br>No Head Support   | Height 21"<br>Shoulder 20"<br>Hip 16" | Height 22"<br>Shoulder 24"<br>Hip 20" | Height 23"<br>Shoulder 28"<br>Hip 24" | Height 24"<br>Shoulder 32"<br>Hip 28" |
| <b>GoLong</b>                       | 42" x 78"                             |                                       | 56" x 78"                             |                                       |
| <b>GoHygiene</b>                    | Chest 32"                             | Chest 36"                             | Chest 40"                             | Chest 44"                             |
| <b>GoComfort</b>                    | Height 23"<br>Shoulder 24"<br>Hip 20" | Height 24"<br>Shoulder 26"<br>Hip 22" | Height 25"<br>Shoulder 28"<br>Hip 24" | Height 26"<br>Shoulder 30"<br>Hip 26" |
| <b>GoActive</b>                     | Chest 40"                             | Chest 44"                             | Chest 48"                             | Chest 52"                             |
| <b>GoReposition</b>                 | 20" x 75"                             |                                       |                                       |                                       |
| <b>GoLimb</b>                       | 15" x 7"                              |                                       |                                       |                                       |
| <b>GoTurn</b>                       | Width 22"                             |                                       |                                       |                                       |
| <b>GoAmputee</b>                    | Height 28"<br>Hip 15"                 | Height 33"<br>Hip 18"                 | Height 38"<br>Hip 22"                 | Height 43"<br>Hip 26"                 |
| <b>GoWalk</b>                       | Hips 15"                              | Hips 18"                              | Hips 22"                              | Hips 26"                              |