

Baby Bassinet

Connor Series



Configuration Examples

2 Door Cabinet



BB-CDR-TM

3 Shelves



BB-CD3-TM

1 Drawer



BB-CD1-TM

No Shelf



BB-CD0-TM

Customize Your Baby Bassinet



Features and Options



Extended Shelf



Flame Retardant Mattress
(26" L x 12" W x 2" H)



Infant Tub
(26.5" L x 13" W x 10" H)

Standard Dimensions

Overall Height (floor to handle)	38.5" (97.8 cm)
Height from Floor to Top of Mattress (floor to baby)	33.25" (84.5 cm)
Overall Width	31" (78.7 cm)
Overall Depth	17.75" (45 cm)
Casters	3" (7.6 cm) 2 locking, 2 swivel