


# Height Adjustable Baby Bassinet

## Jordan Series



**Height Range (floor to mattress) from 28" - 36" (71 cm - 91.4 cm)**



MAXIMUM HEIGHT



LOWEST HEIGHT

# Customize Your Bassinet



## Options



Pull Out  
Writing Surface



Flame Retardant Mattress  
(26" L x 12" W x 2" H)



Infant Tub  
(26.5" L x 13" W x 10" H)



Key  
Lock




Punch Code  
Lock




Electronic  
Lock


## Standard Finishes\*




Hard Rock  
Maple



Candlelight




Dolce  
Vita




Medium  
Oak



Chocolate  
Pear



Snow  
White



Richmond  
Grey

\*Custom colors are available, please contact us for details.

## Available Material

Thermofoil

## Standard Dimensions\*\*

	With Pull Out	Without Pull Out
Overall Height (floor to handle)	29.125" - 37.375" (74 cm - 95 cm)	28.5" - 36.75" (72.4 cm - 93.3 cm)
Height from Floor to Top of Mattress (floor to baby)	28.5" - 36.75" (72.4 cm - 93.3 cm)	27.75" - 36" (70.5 cm - 91.4 cm)
Overall Width	32" (81.3 cm)	32" (81.3 cm)
Overall Depth	21" (53.3 cm)	21" (53.3 cm)
Casters	5" (12.7 cm) 3 locking, 1 swivel	5" (12.7 cm) 3 locking, 1 swivel
Pull Out Work Surface	12.125" w x 13" d (30.8 cm x 33 cm)	n/a

\*\*Custom dimensions are available, please contact us for details.