


Height Adjustable Baby Bassinet

Lily Series



Height Range (floor to mattress) from 26.75" - 35.25" (70 cm - 89.5 cm)



MAXIMUM HEIGHT








LOWEST HEIGHT




Customize Your Bassinet



Options

				
Pull Out Writing Surface	Flame Retardant Mattress (26" L x 12" W x 2" H)	Key Lock	Punch Code Lock	Electronic Lock

Standard Finishes*

						
Hard Rock Maple	Candlelight	Dolce Vita	Medium Oak	Chocolate Pear	Snow White	Richmond Grey

*Custom colors are available, please contact us for details.

Available Material

Thermofoil

Standard Dimensions**

	With Pull Out	Without Pull Out
Overall Height (floor to handle)	30.25" - 36.5" (76.84 cm - 92.71 cm)	29.125" - 37.625" (74 cm - 95.5 cm)
Height from Floor to Top of Mattress (floor to baby)	27.875" - 34.625" (70.8 cm - 87.95 cm)	26.75" - 35.25" (68 cm - 89.5 cm)
Overall Width	37.7" (95.8 cm)	37.7" (95.8 cm)
Overall Depth	21" (53.3 cm)	21" (53.3 cm)
Casters	5" (12.7 cm) 3 locking, 1 swivel	5" (12.7 cm) 3 locking, 1 swivel
Pull Out Work Surface	12.125" w x 13" d (30.8 cm x 33 cm)	n/a

**Custom dimensions are available, please contact us for details.