


Overbed Tables



Various top styles available



Vanity drawer with mirror

Seamless cup holders



Easy, one-hand activated secondary surface



Integrated accessory rail (urinal hook also available)

Shock-absorbent gas spring cylinder for easy height adjustment

OT-F-TYYS-XU

Various base styles available



U Base



H Base



C Base

Lowest base on the market
Exposed steel; no plastic cover allows for easy cleaning

Configuration Examples

Single Top



OT-S-TYYV-XU

Dual Top



OT-F-TYYS-XU

Specialty



JELLYBEAN SHAPE



TABLET MOUNT



OT-S-TYYN-XC



OT-F-TYYS-XU



WHITEBOARD




SOLID SURFACE


Customize Your Overbed Tables




Table Top Options



Standard Top Finishes*



Standard Base Finishes*




*Custom finishes are available, please contact us for details.

Accessories




Flip Table Tops**



**Flips are not available on laminate, storage, jellybean, c-base, or solid surface overbed tables.

Storage Table Options



Storage Table Configuration Examples



Standard Dimensions**

Single Table Top	32" L x 17" W (81.3 cm x 43.2 cm)
Flip Table Top	13.5" L x 17" W (34.3 cm x 43.2 cm)
Storage Table Top	43.5" L x 19" W (110.5 cm x 48.3 cm)
Height Range	28" Low - 42" High (71 cm - 106.7 cm)

**Custom dimensions are available, please contact us for details.

	U Base	H Base	C Base
Overall Length	32.5"	34.38"	30.75"
Overall Width	19"	18.25"	21.25"
Internal Width (between the legs)	14"	29.38"	25"
Overall Height from the Floor	Approx. 1.9"	Approx. 1.9"	Approx. 1.9"