

Care Series

Three Function Electrical Bed

AC-200C



Key Features:

- Easily removable head and foot boards
- Instant CPR release lever
- 320 lb (145 kg) weight capacity
- Bed bumpers for added protection
- Three positional functions (HILO, Fowler and Knee)
- Brake mechanism located at foot of bed
- Battery back-up
- Caregiver-controlled "flat surface" function
- Fowler angle indicator
- 6 IV pole holders
- 2 Foley bag holders
- Linak DC actuators
- Patient control on inner siderails

The AC-200C bed provides caregivers the control they need in the recovery environment. The simplistic design and distinctive features bring efficiency to today's medical and surgical environment.

The Amico Care AC-200C bed features an open architecture for surfaces and removable head and foot boards which provides caregivers easy access in the event of an emergency. The head and foot boards feature an exclusive blend of polypropylene for cleanup and durability.

The functions of the bed are electronically controlled and equipped with a back-up battery that allows caregivers to make adjustments to the bed during transport and unexpected power outages.

Technical Specifications



Model Number

Model Number	AC-200C
--------------	---------

Overall Dimensions

Overall Length	88" (224 cm)
Overall Width	42" (107 cm)
Patient Surface	35" x 77" (89 cm x 196 cm)
Height Range (to the surface of the litter)	
High	29" (74 cm)
Low	18" (46 cm)

Litter Positioning

Backrest	0° - 65°
Knee Gatch	0° - 30°

Weight Capacity & Caster Diameter

Weight Capacity	320 lbs (145 kg)
Caster Diameter	5" (13 cm)

Electronics

Voltage	110 V @ 60 Hz or 220 V @ 50 Hz
---------	--------------------------------

Optional Features & Accessories:

- Chart holder
- 6" (15 cm) foam mattress
- Provision for trapeze
- Handset control

Amico Beds Corporation reserves the right to change product specifications without notice.

Warranty: One year parts only.

NOTE: These products are for Teaching and Training Simulations only.