

Teaching & Training Solutions

Surface Mount Headwalls



Sapphire/Regal Series

Product Benefits:

Amico has developed a revolutionary headwall that is self-contained and alleviates the headaches associated with a nursing simulation lab installation.

The traditional method of outfitting a simulation lab with medical gas services can cost thousands of dollars for pumps, compressors and a network of plumbing. This Plug & Play training headwall simulates Oxygen, Air and even create suction simultaneously via the vacuum outlet by connecting the headwall to a portable compressor.

The compressor can be located anywhere you wish, you can even provide your own in-wall source if required and hook to a Central Compressor with or without a Central Vacuum Pump. The electrical receptacles can have the option of "hard wire".

Product Applications:

Educational Facilities

Teaching/Training

Universities

Simulation Centers

Vocational Schools

Nursing Labs


Mock Up Rooms

Institutions


Standard Colors:



Wild Cherry



Fusion Maple



Antique White




Key Features:

- Plug & Play headwall
- Typical widths: 24" (60.96 cm), 36" (91.44 cm), 48" (121.92 cm), 60" (152.4 cm)
- One full length, scratch-resistant accessory rail
- Three standard colors of hospital grade laminate to choose from
- Super silent, HP air compressor to simulate medical gases
- Rear connection option available for electrical and/or compressor and vacuum
- Mobile floor stand option
- Optional light for units 48" or wider


Technical Specifications




Sapphire Teaching & Training Headwall



Regal Teaching & Training Flatwall



Sapphire with Light Teaching & Training Headwall



Item #	Description	Item #	Description
1	Duplex Receptacle	5	Nurse Call Cover
2	Vacuum Switch	6	Medical Air Outlet
3	Vacuum Outlet	7	Compressed Air Outlet
4	Oxygen Outlet	8	Light (Optional)

* Not Shown

Distributed By:

Amico Beds Corporation reserves the right to change product specifications without notice.

Warranty: One year parts only.