

Teaching & Training Solutions

Nurse Call Simulation Unit



Key Features:

- Functions wirelessly over WiFi
- Nurse Call communication
- Code Blue emergency calls
- Push button controls for standard calls, Code Blue and cancelling a call
- 1/4" Jack for Call cord compatibility
- Easy installation



The Amico Nurse Call Simulation Unit can connect with and make calls wirelessly to a WiFi connected computer or laptop, simulating a Nurse Call station. The intuitive and simplistic design makes it a perfect fit for simulation and classroom environments. The unit can make standard calls as well as Code Blue emergency calls with a simple press of a button. The ease of installation and ease of use makes it a great option for a variety of situations. The Amico Nurse Call Simulation Unit is ideal for bringing the functionality of a Nurse Call station straight to a simulation environment.

Technical Specifications



Model Number



Model Number	N2P-NURSECALL
--------------	---------------

Overall Dimensions

Width	6.50" (165 mm)
Height	4.56" (116 mm)
Depth	2" (51 mm)

Compatibility

2.4 GHz WiFi Connection
Windows 10 or later



- | | |
|------------------------|--------------------|
| 1. Standard Nurse Call | 4. Speaker |
| 2. Cancel Call | 5. Call Cord Jack |
| 3. Code Blue Call | 6. Indicator Light |

Distributed By:

Amico Beds Corporation reserves the right to change product specifications without notice.

Amico Beds Corporation | 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada | 600 Prime Place, Hauppauge, NY 11788, USA
Toll Free Telephone: 1.877.462.6426 | Toll Free Fax: 1.866.440.4986 | Telephone: 905.764.0800 | Fax: 905.764.0862 | www.amico.com

AB-SP-TTS-N2P-NURSE-CALL 10.16.2025