

Teaching & Training Solutions

Universal Receptacle

Sapphire Series



Product Benefits:

Amico has developed a revolutionary headwall that is self-contained and alleviates the headaches associated with a nursing simulation lab installation.

The traditional method of outfitting a simulation lab with medical gas services can cost thousands of dollars for pumps, compressors and a network of plumbing. This Plug & Play training headwall simulates Oxygen, Air and even create suction simultaneously via the vacuum outlet by connecting the headwall to a portable compressor.

The compressor can be located anywhere you wish, you can even provide your own in-wall source if required and hook to a Central Compressor with or without a Central Vacuum Pump. The electrical receptacles can have the option of "hard wire".

Product Applications:

- Educational Facilities
- Teaching/Training
- Universities
- Simulation Centres
- Vocational Schools
- Nursing Labs
- Mock Up Rooms
- Institutions


Standard Colors:




Wild Cherry



Fusion Maple



Antique White



Key Features:

- One full length scratch resistant accessory rail
- Removable panels for access to electrical services
- Accommodates almost any configuration of medical gas outlets, electrical and communication services
- Three standard colors of hospital grade laminate to choose from
- Outlets allow ease of access for annual maintenance
- Typical widths

Technical Specifications

Model Number

Teaching Training Sapphire – Universal Receptacle

BT-HWXX-01T-0-SGEN-UR

Overall Dimensions


With Lights / Without Lights

48"/60"/72"


Technical Drawings



Typical Drawing 48"



Typical Drawing 60"



Amico Beds Corporation reserves the right to change product specifications without notice.

Warranty: One year parts

Amico Beds Corporation | 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada | 600 Prime Place, Hauppauge, NY 11788, USA
Toll Free Telephone: 1.877.462.6426 | Toll Free Fax: 1.866.440.4986 | Telephone: 905.764.0800 | Fax: 905.764.0862 | www.amico.com

TTS-SP-SAPPHIRE-UNIVERSAL-RECEPTACLE 07.05.2024