

# Teaching & Training Solutions

## Surface Mount Panel Wall



### Product Applications:

- Educational Facilities
- Teaching/Training
- Universities
- Simulation Centers
- Vocational Schools
- Nursing Labs
- Mock Up Rooms
- Institutions

### Key Features:

- Simulates Medical Air, Oxygen and Vacuum with one single air compressor source
- Customizable unit height
- Vertical accessory rails




### Product Benefits:

Amico has developed a revolutionary headwall that is self-contained and alleviates the headaches associated with a nursing simulation lab installation.

The traditional method of outfitting a simulation lab with medical gas services can cost thousands of dollars for pumps, compressors and a network of plumbing. This Plug & Play training headwall simulates Oxygen, Air and even create suction simultaneously via the vacuum outlet by connecting the headwall to a portable compressor.

The compressor can be located anywhere you wish, you can even provide your own in-wall source if required and hook to a Central Compressor with or without a Central Vacuum Pump.



Distributed By:

Amico Beds Corporation reserves the right to change product specifications without notice.